

## Kit Contents:

| Cat. No.<br>/ preps          | FABGK000-Maxi<br>(2 preps) | FABGK003<br>(10 preps) | FABGK003-1<br>(24 preps) |
|------------------------------|----------------------------|------------------------|--------------------------|
| Proteinase K powder+         | 11 mg x 2                  | 11 mg x 10             | 11 mg x 24               |
| FABG Buffer                  | 22 ml                      | 110 ml                 | 265 ml                   |
| W1 Buffer* (concentrated)    | 6.5 ml                     | 33 ml                  | 88 ml                    |
| Wash Buffer** (concentrated) | 3 ml                       | 20 ml                  | 40 ml                    |
| Elution Buffer               | 6 ml                       | 30 ml                  | 60 ml                    |
| FABG Maxi Column             | 2 pcs                      | 10 pcs                 | 24 pcs                   |
| Elution Tube (50 ml tube)    | 2 pcs                      | 10 pcs                 | 24 pcs                   |
| User Manual                  | 1                          | 1                      | 1                        |

## Preparation of ProteinaseK solution (20mg/ml) and Wash Buffer for first use:


| Cat. No:                                     | FABGK000-Maxi<br>(2 preps) | FABGK003<br>(10 preps) | FABGK003-1<br>(24 preps) |
|--|----------------------------|------------------------|--------------------------|
| + ddH <sub>2</sub> O volume for Proteinase K | 0.5 ml                     |                        |                          |
| * ethanol volume for Wash Buffer W1          | 2.5 ml                     | 12 ml                  | 32 ml                    |
| **ethanol volume for Wash Buffer W2          | 12 ml                      | 80 ml                  | 160 ml                   |

## Specification:

Principle: spin column (silica membrane)  
 Sample Size : up to 10 ml of fresh/ frozen blood;  
 up to  $1 \times 10^8$  of cultured cells  
 Column Capacity: 500 µg of DNA  
 Average DNA yield : 35 µg/ 1 ml whole blood  
 Handling Time: 1 hour  
 Elution Volume: 0.75~1.5 ml

## Required material to be provided by user

Pipettors and pipet tip  
 Centrifuge: should be capable of producing a force of 4,000 x g  
 Thermal incubator  
 Oven (optional)  
 Ethanol (96~100%)  
 Vortex



## Important Notes:

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
2. Preheat a thermal incubator to 60 °C before the operation.
3. Use a centrifuge with a swinging bucket rotor and a force of 4,000 ~ 6,000 x g for in all centrifugation steps.
4. Add 500 µl of sterile ddH<sub>2</sub>O to proteinase K tube to make a 22 mg/ ml stock solution. Vortex and make sure that Proteinase K powder has been completely dissolved. Store the stock solution at 4 °C.
5. Preheat the Elution Buffer or ddH<sub>2</sub>O for step 11 (Elution step).

## Protocol: for Blood DNA Extraction

### Please Read Important Notes Before Starting The Following Steps.

1. Transfer up to 10 ml sample (whole blood, buffy coat) to a 50 ml centrifuge tube (not provided).  
 --- If the sample volume is less than 10 ml, add PBS to make volume to 10 ml.
2. Add 500 µl of Proteinase K (20 mg/ml) to the sample and mix well by vortexing. Add 10 ml of FABG Buffer to the sample mixture. **Mix thoroughly by pulse-vortexing.**  
 --- Do not add Proteinase K directly to FABG Buffer.
3. Incubate the sample mixture at 60 °C for 15 min to lyse the sample. During incubation, invert the tube every 3-5 minutes.
4. **(Optional):** If RNA-free genomic DNA is required, add 80 µl of 100 mg/ml RNase A (not provided) to the sample mixture and incubate at room temperature for 10 minutes.
5. Add 10 ml of ethanol (96 ~100%) to the sample mixture. **Mix thoroughly by vortexing. If precipitate appears, break it by pipetting.**
6. Place a FABG Maxi Column to a 50 ml centrifuge tube (not provided). And transfer 15 ml of sample mixture (ethanol added) (including any precipitate) carefully to the FABG Maxi Column. Close the cap and **centrifuge at 4,000 ~ 6,000 x g for 3 min.**

7. Discard the flow-through and transfer the rest sample mixture to the same FABG Maxi Column. Close the cap and **centrifuge at 4,000 ~ 6,000 x g for 3 min** and discard the flow-through.
8. Add 4 ml of W1 Buffer (ethanol added) to the FABG Maxi Column. Close the cap and **centrifuge at 4,000 ~ 6,000 x g for 3 min**. Discard the flow-through and place the FABG Maxi Column back in the 50 ml centrifuge tube.
  - Make sure that ethanol has been added into W1 Buffer when first open.
9. Add 7 ml of Wash Buffer (ethanol added) to the FABG Maxi Column. Close the cap and **centrifuge at 4,000 ~ 6,000 x g for 15 min**. Discard the flow-through and place the FABG Maxi Column back in the 50 ml centrifuge tube.
  - Make sure that ethanol has been added into Wash Buffer when first open.
  - **Important Step!** Make sure the residual liquid will be removed completely after centrifugation. It might be necessary to do a further drying by placing the column in a vacuum oven at 70 °C for 3 minutes.
10. Place the FABG Maxi Column into a new 50 ml centrifuge tube. (Elution Tube, provided)
11. Add 0.75 ~1.5 ml of preheat Elution Buffer or ddH<sub>2</sub>O (pH 7.5- 9.0) to the membrane center of the FABG Maxi Column. **Stand the FABG Maxi Column for 5 min at room temperature.**
  - **Important Step!** For effective elution, stand the FABG Maxi Column for 5 minutes is required to make sure Elution Buffer is absorbed completely by column membrane.
  - Standard volume for elution is 0.75 ~1.5 ml. If higher DNA yield is required, repeat the DNA Elution step (step 11) to increase DNA recovery.
12. **Centrifuge at 4,000 x g for 2 minutes** to elute total DNA.

### Protocol: (for Cultured Cell DNA Extraction)

1. Transfer up to  $1 \times 10^8$  of cells to a 50 ml centrifuge tube (not provided). **Centrifuge at 4,000 ~ 6,000 x g for 5 minutes** to pellet the cells. (If using adherent cells, trypsinize the cells before harvesting.)
2. Resuspend the cells with 10 ml of PBS.
3. Follow the Blood protocol starting from step 4.

### Troubleshooting

| Possible reasons   | Solutions  | Possible reasons                                       | Solutions   |
|--|--|--|---|
| <b>Low or no yield of genomic DNA</b>  |  | <b>Elution of genomic DNA is not efficient</b>         |   |
| <b>Poor cell lysis</b>   |  | pH of water (ddH <sub>2</sub> O) for elution is acidic | Make sure the pH of ddH <sub>2</sub> O is between 7.5- 9.0.   |
| Poor cell lysis because of insufficient Proteinase K activity                                | Repeat the extraction procedure with a new sample. Use a fresh or well-stored Proteinase K stock solution.   |  | Use Elution Buffer (provided) for elution.  |
| Poor cell lysis because of insufficient mixing with FABG buffer                              | Repeat the extraction procedure with a new sample. Mix the sample and FABG Buffer immediately and thoroughly by pulse-vortexing.                               |  | Elution Buffer or ddH <sub>2</sub> O is not completely absorbed by column membrane  |
| Poor cell lysis because of insufficient incubation time                                      | Repeat the extraction procedure with a new sample. Extend the incubation time and make sure that no residual particulates remain.                              | <b>Column is clogged</b>                               |   |
| <b>Ethanol is not added into the lysate before transferring into FABG Maxi Column</b>        | Repeat the extraction procedure with a new sample.   | Blood sample contains clots                            | Repeat the extraction procedure with a new sample. Mix the blood sample well with anti-coagulant to prevent formation of blood clots. |
| <b>Incorrect preparation of Wash Buffer</b>  |  | Sample is too viscous                                  | Reduce the sample volume.   |
| Ethanol is not added into W1 and Wash Buffer when first                                      |  | <b>Degradation of eluted DNA</b>                       |   |
| The volume or the percentage of ethanol is not correct before adding into W1 and Wash Buffer | Make sure that the correct volumes of ethanol (96- 100 %) is added into W1 and Wash Buffer when first open. Repeat the extraction procedure with a new sample. | Sample is old  | Always use fresh or well-stored sample for genomic DNA extraction.  |
|  |  | Buffer for gel electrophoresis contaminated with DNase | Use fresh running buffer for gel electrophoresis.   |