


Genaxy Dual Filter Tips



Genaxy dual filter tips are beneficial for setting PCR reactions and handling of radioactive, biotoxic, Chemical corrosive and volatile samples.

Dual layered filter consists of two layers:
Upper filter: upper layer of filter is blue in color having pore size of 5µm. It prevents bio-molecules efficiently
Lower filter: lower layer of filter is white in color having pore size of 10µm. It blocks droplets, splashes and suspended particles.

Unique Features

- » Dual barrier filters with different pore size
- » DNase, RNase, endotoxin and pyrogen free
- » Superb hydrophobicity of polyethylene (UHMW)
- » Maximum air permeability
- » Reinforced blocking property

Cat. Number	Description	Qty. / Pack
GEN-DFB-10	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-10-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-10-R-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-10XL	0.5-20µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-10XL-S	0.5-20µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-10XL-R-S	0.5-20µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-200	1-200µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-200-S	1-200µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-200-R-S	1-200µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-1000	100-1000µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-1000-S	100-1000µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-1000-R-S	100-1000µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	6 racks

Genaxy Scientific Pvt. Ltd.