

## Blotters

### Western Blotters

Catalog No.	Description	Price (Rs.)
SB10	omniPAGE Blot Mini, 10 x 10cm System including blotting insert, tank and lid, 4 cassettes, 16 fibre pads, cooling pack	65,000.00
SB10W	omniPAGE Blot Mini Wide, 20 x 10cm System including tank and lid, 4 cassettes, 16 fibre pads, cooling pack	88,800.00
SB20	omniPAGE Blot MaxiWAVE, 20x 20cm System including tank and lid, 4 cassettes, 16 fibre pads, cooling coil	1,29,600.00
<b>Accessories</b>		
VS10BI	1x Mini Blotting Insert, 4x compression cassettes and 16x fibre pads	39,600.00
VS10BI-HI	1x Mini High Intensity Blotting Insert, 2x compression cassettes and 8x fibre pads	54,270.00
VS10WBI	1x Mini WIDE Blotting Insert, 4x compression cassettes and 16x fibre pads	48,720.00
VS20BI	WAVE Maxi Platinum Wire Blotting Insert, 4x compression cassettes and 18x fibre pads – standard”	63,624.00
VS20BI-HI	WAVE Maxi Plate Electrode Blotting Insert, 1x compression cassette and 6x fibre pads – highintensity”	82,800.00
SB10WC	omniPAGE Blot Mini Wide Cassette	6,000.00
SB20C	omniPAGE Blot Maxi Cassette	7,656.00
SB10F	Fibre Pads – pk/8	4,590.00
SB10WF	Fibre Pads – pk/8	6,000.00
SB20F	Fibre Pads – pk/8	7,660.00



VS10BI

### High Intensity Transfer System

The Sub Blot Mini and WAVE high intensity transfer systems combine the cooling capacity of wet submarine blotting units with the speed of semi-dry transfer systems. Both systems utilize plate electrodes to create a higher strength electric field and greater current density than conventional wire-electrode systems.

Catalog No.	Description	Price (Rs.)
SB10-HI	omniPAGE Blot Mini, 10 x 10cm, High Intensity System including tank and lid, 2 cassettes, 8 fibre pads, cooling pack	90,000.00
SB20-HI	WAVE Sub Blot, 20 x 20cm, High Intensity System including tank and lid, 1 cassettes, 6 fibre pads, cooling coil	1,32,000.00
<b>Accessories</b>		
VS10BI-HI	omniPAGE Blot Mini Insert – includes 2 cassettes and 8 fibre pads	54,270.00
SB10C	Sub Blot Mini Cassette	4,590.00
SB10F	Fibre pads - Pack of 8	4,590.00
VS20BI-HI	WAVE Maxi Plate Electrode Blotting Insert, 1x compression cassette and 6x fibre pads – highintensity	82,800.00
SB20C	Sub Blot Maxi Cassette	7,660.00
SB20F	Fibre pads - Pack of 8	7,660.00




SB20-HI

### Tank Sub Electro Blotters

Tank Sub Electroblotters can utilize up to five gel blot cassettes at any one time for proteins. Accommodates gel thickness from 0.25 up to 3 mm and is especially useful in high throughput laboratories.

Catalog No.	Description	Price (Rs.)
EBM10	Electro Blot Mini, 10 x 10cm System for five cassettes, with tank and lid, 5x cassettes, 24x fibre pads and cooling pack	78,000.00
EBM20	Electro Blot Maxi, 20 x 20cm Five Position, including 5 cassettes, with Tank & Lid, 24 fibre pads, cooling pack	1,08,000.00
<b>Accessories</b>		
SB10C	Sub Blot Mini Cassette	4,590.00
SB10F	Fibre pads - Pack of 8	4,590.00
SB20C	Sub Blot Maxi Cassette	7,660.00
SB20F	Fibre pads - Pack of 8	7,660.00



Tank Sub Blotter

Unit dimension (WxHxD)	Mini:19x13x19cm Maxi:24x16x26cm
Max.sample capacity	
Mini	5 Blots, 10x10cm
Maxi	5 Blots, 20x20cm 20 Blots, 10x10cm
Buffer volume	
Mini	Min 1000ml; Max 1500ml
Maxi	Min 4300ml; Max 6000ml



**Semi Dry Blotter**

Max. sample capacity	
Mini	1 Blot, 8x8.5cm
Maxi	1 Blot, 16x17.5cm, 4 Blots, 8x8.5cm, 2 Blots, 16x8.5cm
Maxi Plus	1 Blot, 33x45cm
Maxi Long	1 Blot, 20x50cm
Buffer volume	Mini:5ml-Maxi:20ml, Maxi Plus:75 ml-Maxi Long:50 ml



**Dot & Slot Blotter**



**MicroDOC**

## ☉ Semi Dry Blotters

These Semi Dry Blotters offer rapid transfer times for DNA, RNA and protein blotting – typically 15 to 30 minutes. All units can be used for all types of blotting: western, Southern and northern blotting. Screw down lid- gels from 0.25 up to 10mm thick can be blotted.

Catalog No.	Description	Price (Rs.)
SD10	Semi Dry Mini, 10 x 10cm System	69,600.00
SD20	Semi Dry Maxi, 20 x 20cm System	1,44,000.00
CSL-SDCAB	Semi dry Blotter Cables (Red & Black)	6,240.00

## ☉ Dot & Slot Blotters

Four different sample number size and style of Hybridisation Manifold are offered -two types of Dot Blotter and two types of Slot Blotter. Typical applications include clone screening with DNA / RNA probes in Southern/northern blots and immunological screening with antibodies in western blots..

Catalog No.	Description	Price (Rs.)
CSL-D48	Dot Blotter, 48 Sample	63,600.00
CSL-D96	Dot Blotter, 96 Sample	73,200.00
CSL-S24	Slot Blotter, 24 Sample	73,200.00
CSL-S48	Slot Blotter, 48 Sample	78,000.00

## Gel Documentation Systems

### ☉ microDOC™

MicroDOC is the researcher's choice for a simple ultra-compact gel documentation system with Computer-free operation. Supplied with a 4GB storage card and 55mm ethidium bromide filter as standard, while an optional SYBR filter is also available. Having 18 Megapixel digital camera with improved high sensitivity CCD sensor and processor. Image visualized within a large 8" TFT colour monitor. Light weight compact hood with easy access door and built-in inner lights..

Catalog No.	Description	Price (Rs.)
CSL-MICRODOC	Compact Gel documentation system without printer and without transilluminators	4,44,000.00
CSL-MDOCUV312	microDOC System with UV Transilluminator (UVT312)	POR
CSL-MDOCUV254	microDOC System with UV Transilluminator (UVT254)	POR
CSL-MDOCUV365	microDOC System with UV Transilluminator (UVT365)	POR
CSL-MDOCUV254/312	microDOC System with UV Transilluminator (UVT254/312)	POR
CSL-MDOCUV254/365	microDOC System with UV Transilluminator (UVT254/365)	POR
CSL-MDOCUV312/365	microDOC System with UV Transilluminator (UVT312/365)	POR
CSL-MDOCUV3121D	microDOC System with UV Transilluminator (UVT312) & 1D Software	POR
CSL-MDOCUV2541D	microDOC System with UV Transilluminator (UVT254) & 1D Software	POR
CSL-MDOCUV3651D	microDOC System with UV Transilluminator (UVT365) & 1D Software	POR
CSL-MDOCUV254/3121D	microDOC System with UV Transilluminator (UVT254/312) & 1D Software	POR
CSL-MDOCUV254/3651D	microDOC System with UV Transilluminator (UVT254/365) & 1D Software	POR
CSL-MDOCUV312/3651D	microDOC System with UV Transilluminator (UVT312/365) & 1D Software	POR

#### Accessories

CSL-PRINT	Thermal Printer	2,47,200.00
CSL-MDOCEB	Microdoc ethidium bromide filter	74,400.00
CSL-MDOCSBRG	Microdoc SYBR filter	74,400.00
CSL-PRTPAP	Replacement printer paper	6,720.00
CSL-UVSCRN	UV to white light conversion screen for transilluminator	73,720.00
CSL-MDOCWLB	White light box for Micro doc	57,600.00