

Gel Cutting Tips

GeIX Disposable Gel Excision Tips

GeIX disposable gel excision tips offer a unique and convenient way to excise bands from gels. GeIX tips allow completely safe and contamination free removal of the gel band of interest without the need to handle potentially dangerous scalpels. The gel band is neatly extracted and contained in the tip, when it is then dispensed into a tube for further analysis. Once used, the tip is ejected for disposal.

- + Eliminates scalpel damage to transilluminator or gel tray.
- + Designed for removing bands from agarose and protein gels.
- + Push button gel and tip release mechanism.
- + Fits standard 1000 μ l pipettors.
- + DNase and RNase free & Autoclavable. Sterile versions are also available

Catalog No.	Description	Price (Rs.)
CSL-GELX4	4mm x 1mm, Gel Cutting Tips, 250 Tips/ Bag	14,400.00
CSL-GELX4RACK	4mm x 1mm, Gel Cutting Tips, 5 Racks of 48 Tips	28,800.00
CSL-GELX6.5	6.5mm x 1mm, Gel Cutting Tips, 250 Tips/ Bag	14,400.00
CSL-GELX6.5RACK	6.5mm x 1mm, Gel Cutting Tips, 5 Racks of 48 Tips	28,800.00

Special Applications Electrophoresis

Clinical Electrophoresis/Cellulose Acetate Electrophoresis

- + The CellasGEL range includes both equipment and consumables to assist in the research and diagnosis of specific disease states.
- + The ideal tank for standard 'dry' membrane and 'wet' gel cellulose acetate techniques. Two adjustable supports, which can be positioned anywhere within the tank, readily accommodate different lengths of dry cellulose acetate membrane to a maximum 20cm

Catalog No.	Description	Price (Rs.)
CSLCELLAS	Horizontal Unit for Cellulose Acetate Electrophoresis	81,600.00

Recommended Accessories

CSLDENS	TurboScan Software Densitometer (Excludes Computer and Scanner)	POR
CSLSCAN	Flatbed scanner for TurboScan software	1,08,000.00
CSLGEL2.514250	2.5x14 CellasGEL 250 μ m, standard thickness, 100/pack	17,400.00
CSLBDG8.5S	1x8.5cm bridge for 1x5.7x14cm or 2x2.5x14cm CellasGEL strips	11,358.00
CSLAPPS22	1x2-specimen semi-micro applicator for 2.5x14cm Strip	61,200.00

Cellulose Acetate Strips

Catalog No.	Description	Price (Rs.)
CSLGEL2.514200	2.5x14cm, CellasGEL 200 μ m, standard thickness, 100/Pack	15,916.00
CSLGEL2.514190	2.5x14cm, CellasGEL 190 μ m, high resolution, 100/Pack	17,400.00
CSLGEL2.517200	2.5x17cm, CellasGEL 200 μ m, 25/Pack	9,000.00
CSLGEL5.714500	5.7x14cm, CellasGEL 500 μ m, high volume, 25/Pack	10,200.00
CSLGEL5.714250	5.7x14cm, CellasGEL 250 μ m, standard thickness, 25/Pack	9,000.00
CSLGEL5.714200	5.7x14cm, CellasGEL 200 μ m, standard thickness, 25/Pack	9,000.00
CSLGEL5.714190	5.7x14cm, CellasGEL 190 μ m, high resolution, 25/Pack	9,000.00
CSLGEL1717200	17x17cm, CellasGEL 200 μ m	19,440.00

Sample Applicators

Catalog No.	Description	Price (Rs.)
CSLAPPS4SP	1x 4-specimen semi-micro applicator for 5.7x14cm	61,200.00
CSLAPPS4	1x 4-specimen semi-micro applicator for 5.7x14cm	61,200.00
CSLAPPS6	1x 6-specimen semi-micro applicator for 5.7x14cm	61,200.00
CSLAPP8	1x 8-sample micro applicator for 5.7x14cm	69,600.00



Cellulose Acetate Electrophoresis



Bridges



Applicator

Kits


Catalog No.	Description	Price (Rs.)
CSLKITCU	CellasKIT: 1-patient test for serum and concentrated urine IFE 2x semi-micro tests (1x serum and 1x urine) for 1 patient performed on 6xCellasGEL strips (2.5x14cm; 190 μ m) on 3x bridges in 1x CSL-CELLAS. Kit content, sufficient for 5 patients (10x semi-micro tests): 30 CellasGEL strips & TGS buffer; Coomassie stain, clearing & saline solutions; volumetric distributors & antisera (anti- IgG, IgA, IgM, Igk & Ig λ); blotting paper & mylar film. Excludes: Destain.	76,800.00
CSLKITI2432	CellasKIT: 6-patient test for serum IFE 6x semi-micro or 8x micro tests (1/patient) on 6x CellasGEL strips (5.7x14 cm;190 μ m), performed on 6x bridges in 2x CSL-CELLAS units. Kit content, sufficient for 24 patients (semi-micro) & 32 (micro): 24 CellasGEL strips & Tris- Hippurate buffer; Amidoblack stain, saline & clearing solutions; volumetric distributors & antisera (anti- IgG, IgA, IgM, Igk & Ig λ); blotting paper & mylar film. Excludes: Destain.	92,400.00
CSLKITSP100200	CellasKIT: serum proteins 4x semi-micro or 8x micro tests for each CellasGEL strip (5.7x14cm; 200 μ m); 3 strips / CSL-CELLAS performed on 3x bridges. Kit content, sufficient for 100x semi-micro or 200x micro tests: 25 CellasGEL strips & Tris-Hippurate buffer; Ponceau S stain, destain & clearing solutions; blotting paper & mylar film.	30,600.00
CSLKITSP150200	CellasKIT: serum proteins (high resolution) 6x semi-micro or 8x micro tests for each CellasGEL strip (5.7x14cm; 190 μ m); 3 strips / CSL-CELLAS performed on 3x bridges. Kit content, sufficient for 150x semi-micro or 200x micro tests: 25 CellasGEL strips & Tris-Hippurate buffer; Coomassie stain, citric acid & clearing solutions; blotting paper & mylar film.	46,800.00
CSLKITHG100	CellasKIT: haemoglobin 4x semi-micro tests for each Cellas GEL strip (5.7x14cm; 250 μ m); 3 strips / CSL-CELLAS performed on 3x bridges. Kit content, sufficient for 100x semi-micro tests: 25 CellasGEL strips & Tris-Glycine buffer; Ponceau S stain, destain & clearing solutions; blotting paper & mylar film."	34,800.00
CSLKITLP100	CellasKIT: lipoproteins 4x semi-micro tests for each CellasGEL strip (5.7x14cm; 250 μ m); 3 strips / CSL-CELLAS performed on 3x bridges. Kit content, sufficient for 100x semi-micro tests: 25 Cellogel strips & Tris- Hippurate buffer; Sudan Black stain & clearing solution; blotting paper, mylar film and 1 plastic box.	34,800.00

*Dry Membranes & Applicators are also available.


Comet Assay Tanks

- + Cleaver Scientific COMET assay tanks are specialized units to study single cell gel electrophoresis (SCGE), a technique made popular by drug toxicology and carcinogenesis studies for the detection and quantitation of DNA damage in cells.
- + Each tank's robust construction from ebony acrylic ensures that cells remain free of exposure to background light and DNA damage during electrophoresis, while a cooled central platform provides a convenient surface for slide preparation and control of slide temperature during the assay.
- + Available in three slide formats.

Catalog No.	Description	Price (Rs.)
CSL-COM10	Comet Assay Tank for 10 Slides	56,982.00
CSL-COM20	Comet Assay Tank for 20 Slides	66,024.00
CSL-COM40	Comet Assay Tank for 40 Slides	92,400.00



Celas Boxes



Comet Assay Tank