


DNA Clean-Up Kits

FavorPrep™ Gel Purification Kits


- ❑ Rapid extraction of DNA fragments from agarose gels.
- ❑ Highly pure DNA suitable for even PCR and sequencing.
- ❑ Sample size : Upto 200mg agarose gel.
- ❑ No Phenol chloroform extraction, CsCl gradients or ethanol precipitation.
- ❑ 70 bp to 12kb DNA fragment can be purified. The recovery rate is upto 85%.
- ❑ It take about 30 minutes.



CATALOG NO.	Description	Kit Size
FAGPK 001	FavorPrep™ Gel Purification Kit	50 Preps
FAGPK 001 -1	FavorPrep™ Gel Purification Kit	200 Preps
FAGPK 001-2	FavorPrep™ GEL Purification Mini Kit	300 Preps
FAMGK 001	FavorPrep™ MicroElute GEL Purification Kit	50 Preps
FAMGK 001-1	FavorPrep™ MicroElute GEL Purification Kit	200 Preps
FAGKE 96004	FavorPrep™ 96-well GEL Purification Kit	96 x 4 Plates

FavorPrep™ PCR Clean Up Kits

- ❑ Rapid extraction of DNA fragments from PCR solution.
- ❑ Just 3 easy and quick steps to isolate pure PCR product.
- ❑ Sample size : Upto 100µl PCR Product.
- ❑ No Phenol chloroform extraction, CsCl gradients or ethanol precipitation.
- ❑ 100bp to 12kb DNA fragment can be purified. The recovery rate is upto 95%.
- ❑ It take about 15 minutes



Catalog No.	Description	Kit Size
FAPCK 001	FavorPrep™ PCR Clean-up Kit	50 Preps
EAPCK 001 -1	FavorPrep™ PCR Clean-up Kit	100 Preps
FAPCK 001-2	FavorPrep™ PCR Clean-Up Mini Kit	300 Preps
FAMPK 001	FavorPrep™ MicroElute PCR Clean Up Kit	10 Preps
FAMPK 001-1	FavorPrep™ MicroElute PCR Clean Up Kit	50 Preps
FACKE 96004	FavorPrep™ 96-well PCR Clean-Up Kit	96 x 4 Plates