

# Dry Bath Incubator

## GENAXY Mini Dry Bath Incubator

The MiniT series dry bath incubator combine a microprocessor-controlled, a fast heating micro-tube block. The MiniT series are suitable for a range of uses in life science, molecular biology, environmental and industrial laboratories and a variety of applications including incubation, enzyme reactions, immuno-assays. The small form factor allows these to find a place on even the most crowded laboratory benchtop and easy to move.

### Features

- ❖ LCD display. It is easy to set up and use
- ❖ Small footprint: MiniT series are a personal sized that fits in the palm of your hand
- ❖ Simultaneous display of set and actual time, temperature
- ❖ Temperature deflection adjusting meet with different demands
- ❖ Beep-signal/Stop after program completed
- ❖ Interchangeable aluminum alloy sample blocks provide versatility and allow easy cleaning and disinfecting
- ❖ Mini-dimensions : Mini Dry Bath is a personal sized that fits in the palm of your hand
- ❖ Many aluminum blocks are available covering a range of tube sizes from 0.2ml to 50ml and other blocks



GEN-MiniT



GEN-MiniT-100



GEN-MiniT-100H

### Technical Specification

Model	GEN-MiniT	GEN-MiniT-100	GEN-MiniT-100H
Temperature Control Range	RT+5°C~80°C	RT+5°C~100°C	RT+5°C~100°C
Timing Range	1~999min or 1~999sec	1~999min or 1~999sec	1~999min or 1~999sec
Temp. Control Accuracy	±0.5°C	±0.5°C	±0.5°C
Display Accuracy	±0.1°C	±0.1°C	±0.1°C
Heating Time(25°C to Max.Temp.)	≤12min	≤15min	≤15min
Cooling Time(25°C to 0°C)	----	----	----
Temp. Control Accuracy @40°C	±0.3°C	±0.3°C	±0.3°C
Heating Lid	----	----	Heating Lid
Heating Lid Temperature Range	----	----	RT.+5°C~100°C
Dimension(L x W x H) mm	150 x 110 x 100	152 x 120 x 112	152 x 120 x 112
Net Weight	0.5kg	0.8kg	0.9kg
Power	DC12V	DC24V	DC24V



### Accessories

Catalog No	Type	Description
AS-01131-01	A-MiniT	40x0.2ml
AS-01131-02	B-MiniT	24x0.5ml
AS-01131-03	C-MiniT	15x1.5ml
AS-01131-04	D-MiniT	15x2.0ml
AS-01131-05	E-MiniT	8--12.5x12.5x32mm
AS-01131-06	F-MiniT	4x15ml
AS-01131-07	G-MiniT	2x50ml

### Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot No. 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com