

B.O.D INCUBATOR

FEATURES

- ❖ B.O.D. Incubators are suitable for preservation of vaccines, insulin, liver extracts, chemicals etc. and to make biochemical oxygen demand determination.
- ❖ The lower temperature incubator is made of double walled M.S. Sheet finished with epoxy based powder coated paint. Inner chamber is made of highly polished stainless steel/anodized aluminum.
- ❖ Two doors are provided. Inner door is made of transparent acrylic for inspecting the specimens; outer door is insulated and is fitted with magnetic tape with lock and key.



B O D INCUBATORS are suitable for preservation of vaccines, insulin, liver extracts, chemicals etc. and to make biochemical oxygen demand determination.

CONSTRUCTIONS The low temperature incubator is made of double walled M.S.sheet finished with epoxy based powder coated paint, inner chamber is made of highly polished stainless steel. Two door are provided. inner door is made of transparent acrylic for inspecting specimens outer door is insulated and is fitted with magnetic tape with look and key.

Temperature range from 5 c to 50 c with an accuracy of + 1 c provided with two air circulating fans, high performance compressor, cooling coils, heating elements, digital temperature control and display is also fitted a door operated illumination and a caster wheel for easy mobility.

Cat. No.	Volume Cu. Ft.	Height	Width	Depth
1125/05	6.1	825	505	415
1125/10	8	850	600	500
1125/15	10	880	550	550
1125/20	12	900	650	580
1125/25	17			

B.O.D. INCUBATOR provided with 304 grade stainless steel outer body with P.I.D. controller with display

Cat. No.	Volume Cu. Ft.	Height	Width	Depth
1130/05	6.1	825	505	415
1130/10	8	850	600	500
1130/15	10	880	550	550
1130/20	12	900	650	580
1130/25	17			

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com