

DOUBLE WALL WATER BATH

Double walled, inside made of Stainless Steel 304 Qlty Outside Mild Steel sheet painted in Epoxy Powder Coating. The gap of 25mm between the outer and inner walls insulated with the mineral wool and top of the bath and the concentric rings are made of stainless steel 304 quality. The temperature is controlled by a thermostat from few °C above ambient to 100 ° C with a variation of plus minus 0.5 ° C.

The bath is provided with two indicating lamps, concentric rings constant water level arrangement ON/OFF Switch & wire and plug but without thermometer.

OPTION:

- Option - 1: Stirrer with 1/20 H.P. high speed motor (Mandatory)
- Option - 2: Micro-processed Digital Temperature Controller Cum Indicator
- Option - 3: Micro-processed PID controller having 4 ramp & 4 soak profile
- Alarm: it is possible to provide electronic alarm, if temperature deviates more than ± 5 °C to cost extra.
- Option-1 Is compulsory to be attached with Option -2
- Option-1 Is compulsory to be attached with Option-3

(AN ISO 9001: 2008 CERTIFIED& CE MARKED)



CATALOG NO.	SIZE	NO. OF HOLES	CAPACITY
GEN-WB128-1	300 X 250 X 100 mm	6 HOLES OF 75 mm DIA	8 Ltrs.
GEN-WB128-2	355 X 405 X 100 mm	12 HOLES OF 75 mm DIA	15 Ltrs.

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com