

## GENAXY Incubators

### CONSTRUCTION

- Double walled inside aluminum/ Stainless Steel 304 quality.
- Outside of Mild Steel Sheet furnished in Epoxy Powder Coating.
- Insulation by mineral wool of 75mm door has double glass fitted for observation of specimen
- Heating elements are at bottom and sides.

### TEMPERATURE

- Electronic Digital Temperature Controller cum Indicator from few deg.°C above ambient 70 deg. C with +/- 0.5 deg. C
- Perforated shelves adjustable with Air cooling Fan.
- Option; Automatic timer extra cost .
- Work on 220/230 volts

### (AN ISO 9001: 2008 CERTIFIED & CE MARKED)

- Microprocessed PID controller having 4 ramp & 4 soak profile having set mode, alarm parameter, configuration parameter, profile parameter, set point profile, Pv error indicator , parameter for locking etc.
- Alarm : It is possible to provide electronic alarm, if temperature deviates more than  $\pm 5$  °C to cost extra
- Automatic timer electronic 999 mts to cost extra



## Technical Specification

Catalog No.	Inner Chamber	Size	Capacity	Qty.
GEN-INCA-1	Aluminium	300 X 300 X 300 mm	28	each
GEN-INCS-1	Stainless Steel	300 X 300 X 300 mm	28	each
GEN-INCA-2	Aluminium	355 X 355 X 355 mm	45	each
GEN-INCS-2	Stainless Steel	355 X 355 X 355 mm	45	each
GEN-INCA-3	Aluminium	455 X 455 X 455 mm	95	each
GEN-INCS-3	Stainless Steel	455 X 455 X 455 mm	95	each
GEN-INCA-4	Aluminium	455 X 455 X 605 mm	125	each
GEN-INCS-4	Stainless Steel	455 X 455 X 605 mm	125	each
GEN-INCA-5	Aluminium	605 X 605 X 605 mm	224	each
GEN-INCS-5	Stainless Steel	605 X 605 X 605 mm	224	each
GEN-INCA-6	Aluminium	605 X 455 X 910 mm	252	each
GEN-INCS-6	Stainless Steel	605 X 455 X 910 mm	252	each
GEN-INCA-7	Aluminium	605 X 605 X 910 mm	336	each
GEN-INCS-7	Stainless Steel	605 X 605 X 910 mm	336	each

### Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot No. 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com