

MUFFLE FURNACE**GENAXY MUFFLE FURNACE**

- Made of thick Mild Sheet having Powder Coated
- Ceramic Fiber blanket for high density better insulation which makes it light weight. also
- Provided with Stainless Steel 310 Qlty. Base plate.
- Heating Elements are of Kanthal wire
- Temperature : continuously upto 1150° C, 1200° C intermittent
- Temp. controlled by Digital Indicator Cum Controller
- Work on 220 / 440 Volt AC 50 Hz.
- Option: PID Controller having 4 ramp & 4 soak profile can be fitted in above instruments at an extra cost .
- Alarm: It is possible to provide electronic alarm if, temperature deviates more than +/-5 deg.c to cost extra.

Catalog No.	Size	Capacity	Qty.
GEN-MF-1	100 X 100 X 225MM	1.5 KW	each
GEN-MF-2	125 X 125 X 300MM	3.0 KW	each
GEN-MF-3	150 X 150 X 300MM	3.0 KW	each
GEN-MF-4	200 X 200 X 300MM	3.0 KW	each
GEN-MF-5	200 X 200 X 455MM	3.0 KW	each
GEN-MF-6	200 X 200 X 605MM	3.0 KW	each
GEN-MF-7	300 X 300 X 300MM	9.0 KW	each
GEN-MF-8	300 X 300 X 455MM	12.0 KW	each
GEN-MF-9	300 X 300 X 605MM	18.0 KW	each



(AN ISO 9001: 2008 CERTIFIED& CE MARKED)

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com