

GENAXY MX 100-4A Shaker

The MX100-4A is designed for shaking ELISA plate(96 well plates), TCPS and 1.5/2.0ml centrifuge tubes.

Features

- Variable speed from 200rpm to 1,500rpm
- DC brushless motor for gentle and powerful vibration
- Continuous or timing operation
- Orbital of 3mm for efficient mixing in small sample vessels
- Optional rack for 1.5ml centrifuge tube
- Simple spring holders for easy and quick plate installation



TECHNICAL SPECIFICATION

Model	GEN-MX100-4A
Mixing Orbit	3mm
Speed Setting Resolution	10rpm increment
Speed	200~1500 rpm @Standard 200~1000 rpm @Deep well plate or tube holder
Timer	1min~99h59min
Power	220V / 110V, 2A
Dimensions (mm)	290 x 270 x 120
Net Weight	6.5kg

ACCESSORIES

Code	Description
AS-02041-01	20x1.5ml tube holder (apolegamy)
AS-02041-02	20x0.5m tube holder (apolegamy)

Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com