

GENAXY Mini-6K & 10K Centrifuge

This new Mini Centrifuge accommodates 6 x 1.5-2.0ml tubes & 2 x 8 x 0.2ml strip tubes. This unit is ideal for quick spin downs, microfiltration & cell separations available in 6,000 rpm and 10,000 rpm.

- **Two Models to choose from**
- **Minifuge for 1.5-2.0ml tubes and 8 x 0.2ml strip tubes (available in 6K model only)**

Technical Specification

Model	GEN-MINI-6K	GEN-MINI-10K
Power requirement	30W	40W
Relative centrifugal force (RCF)	2,000 x g	5,000 x g
Max. rotational speed	6,000 rpm	10,000 rpm
Rotors	6x1.5/2.0ml angle rotor 2 - 8x0.2ml strip rotor Rotor for slide	6x1.5/2.0ml angle rotor
Dimensions	175 x 148 x 118 mm	175 x 148 x 118 mm
Weight	0.7kg	0.7kg



Accessories

Code	Description
AS-08061-01	Allen wrench
AS-08061-02	1.5/2.0ml rotor
AS-08061-03	Strip tube rotor
AS-08061-04	Slide rotor
AS-08061-05	Adapters for 0.5-0.6ml tubes
AS-08061-06	Adapters for 0.2ml tubes



GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Plot No. 32-57, Community Centre, Behind Janak Cinema,
Janak Puri, New Delhi-110058, Phone:011-4761999 (100 lines), Fax:011-25534163, Email:sales@genaxy.com