

GENAXY Mini Vortex Mixer

The Mini Vortex is an essential unit in all labs which require quick mixing of individual test tubes (of varying sizes).

Features:-

- **Tiny Vortex with the small footprint**
- **Available in fixed and variable speed versions**
- **Brushless DC motor for maintenance free long life**
- **Deliver upto 4500 RPM for thorough mixing & quick vortexing**
- **Orbital Diameter of 4.5mm**
- **Smart and Robust design - press activation with suspension**
- **Robust construction with Die cast Al. base and rubber pads for stable and silent operation**
- **Inbuilt counterbalance for stable operation**
- **Adjustable speed**



Catalogue No.	GEN-iSwix-Mini
Maximum RPM	4500
Orbital Diameter	4.5 mm
Motor Type	Brushless DC
Speed setting	Yes
Timer setting	No
Maximum Load Capacity	100gm
Maximum Tube Diameter	30mm
Press activation	Yes
Continuous mode	Yes
Protection class	IP 43
Dimensions(W x D x H)	100 x 100 x 67 mm
Weight	0.55 Kg

Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com