

Dry Bath Incubator

GEN-MiniT-H2C/GEN-MiniT-H3

GENAXY MiniT-H2C / MiniT-H3 are 3-in-1 mini dry bath incubator. The MiniT-H2C is a heating and cooling incubator and the MiniT-H3 is a heating only incubator with 3 individually controlled mini dry bath. Many aluminum blocks are available covering a range of tube sizes from 0.2ml, 0.5ml, 1.5ml or 2.0ml tubes blocks. These are suitable for labs with heavy work load.

Features

- ✓ Working area dimension: 3 individually controlled
- ✓ LCD display. It is easy to setup and use
- ✓ Simultaneously display the set and actual time, temperature, Auto-stop and sound-alert after program completed
- ✓ MiniT-H3 has a heating function only and its temperature ranges from RT+5~100°C
- ✓ MiniT-H2C has cooling/heating functions and its temperature ranges from 0°C~99.9°C
- ✓ Temperature deflection adjusting meet with different demands



Technical Specification

Model	GEN-MiniT-H2C			GEN-MiniT-H3		
Working area	Block-A	Block-B	Block-C	Block-A	Block-B	Block-C
Temperature Control Range	0°C~99.9°C	RT.+5~99.9°C		RT.+5~99.9°C		
Heating Time(25 to 99.9°C)			≤20min			
Cooling Time(25 to 0°C)	≤25min (*)	----	----	----	----	----
Temperature Control Accuracy			≤ 0.5°C			
Temp. Control Accuracy @37°C			0.2°C			
Display Accuracy			0.1°C			
Timing Range			1~999h or 1~999 min or 1~999sec			
Block Quantity			3			
Power			DC12V,120W			
Dimension(mm)			360 × 180 × 150			
Net Weight			3kg			

Accessories

Catalog No.	Description
AS-01131-01	A-MiniT, 40x0.2ml
AS-01131-02	B-MiniT, 24x0.5ml
AS-01131-03	C-MiniT, 15x1.5ml
AS-01131-04	D-MiniT, 15x2.0ml
AS-01131-05	E-MiniT, 8--12.5 x 12.5 x 32 mm
AS-01131-06	F-MiniT, 4x15ml
AS-01131-07	G-MiniT, 2x50ml



Genaxy Scientific Pvt. Ltd :

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com