

GENAXY Orbital Shaker

GENAXY OS-200/100 Orbital Shaker is a powerful variable speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control. The shaker can be applicable in different laboratories: in microbiology, chemistry, immunology, biochemistry, and molecular biology



GEN-OS-200



GEN-OS-100

Technical Specification

Model	GEN-OS-200	GEN-OS-100
Shaking Speed	RT+5°C~70°C	RT+5°C~70°C
Orbit	20 mm	10 mm
Independent timer	1min~99h59min	1min~99h59min
Load capacity	2.5kg	3.0kg
Motor	DC brushless	DC brushless
Dimension(L x W x H) mm	280 x 270 x 110	280 x 270 x 110

Features

Digital speed control

Universal platform, flat platform or platforms with clamps

Orbital Shaker can be used in a cold room or incubator

Operating ambient: temperature range +5°C to +40°C

Overall dimensions: 280 x 270 x 110 mm

Power supply: DC24V

Weight: 7.2kg

Compatible in cold room usage

Accessories

S No.	Code	Description
1	AS-04021-01	PP-1 Flat platform with non-slip rubber mat
2	AS-04021-02	PP12-100 Platform with clamps for flasks, 100-150ml
3	AS-04021-03	PP06-250 Platform with clamps for flasks, 250-300ml
4	AS-04021-04	UP-1 Universal platform with adjustable bars for different types of flasks, bottles
5	AS-04021-05	MP-1 Platform with springs

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com