

GENAXY HORIZONTAL SYSTEM

MIDI Horizontal Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.



GENECO-H-MIDI

Technical Specification

Model	GENECO-H-MIDI
Principal Material	Acrylic
Dimensions (L X B X H)	27 X 16 X 8 cm
No. of Samples	13
Gel Size	13 x 13 cm
Trays (L x W)	13 X 13 cm (1), 13 X 6 cm (2), 6.5 X 6 cm (4)
Combs	13 Well x 1.5mm (1) Acrylic 8 Well x 1.5mm (4) Acrylic 3 Well Preparative x 3mm (1) Acrylic
Gel Casting Tray	1 Universal
Platinum Electrode	Red & Black, 1 each
Connecting Cord	Red & Black, 1 each
Lid	1 No.

Horizontal System Includes

Basic Unit (1 No.)

Platinum Electrode Assembly Removable (2 Nos.)

3 Well Preparative Comb 3.0 mm (1 No.)

8 well combs 1.5 mm (4 Nos.)

Gel Running Tray 130 X 130 mm (1 No.)

Gel Running Tray 65 X 60 mm (4 Nos.)

Gel Platform 130 X 65 mm (2 Nos)

Lid (1 No.)

Connecting Cord (2 Nos.)

13 Well Combs 1.5 mm (1 No.)

Accessories

S No.	Catalog No.	Description
1	GENECO-H-MIDI-6	Gel Running tray 6 x 13 cm (2/pack)
2	GENECO-H-MIDI-13	Gel Running tray 13 x 13 cm (1/pack)
3	GENECO-H-MIDI-13UV	UV transparent Gel Running tray 13 x 13 cm (1/pack)
4	GENECO-H-MIDI-GS13	Gel Scoop for 13 X 13 cm gel
5	GENECO-H-LT-15	Levelling table for gel casting 15 x 15 cm
6	GENECO-H-MIDI-EL	Spare electrode 2 Nos.

Combs

S No.	Code	Description
1	GENECO-H-MIDI-1.5-8	Acrylic Comb MINI 8 Well 1.5mm
2	GENECO-H-MIDI-1.5-13	Acrylic Comb MINI 13 Well 1.5mm
3	GENECO-H-MIDI-3.0-3	Acrylic Comb MINI 3 Well 3 mm
4	GENECO-H-MIDI-CS	Comb Stand

Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com