

GENAXY HORIZONTAL SYSTEM

MIDIprep Horizontal Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.



GENECO-H-MIDIprep

Technical Specification

Model	GENECO-H-MIDIprep
Dimensions(L X W X H)	370 x 180 x 100 mm
Principal Materials	Acrylic
Max Gel Size	250x130mm
Gel Platforms	250 x 130mm (1 NO)
13 well 1.5mm	Two
3 well 3mm	One
Platinum Electrodes	Positive & Negative
Gel Casting Tray	One
Connecting Cords	Red & Black
Buffer Required	500ml
Power Required	250V, 100mA

Horizontal System Includes

Basic Unit (1 No.)

Platinum Electrode Assembly Removable (2 Nos.)

3 Well Preparative Combs 3.0 mm

3 Well Preparative Combs 6.0 mm (1 No.)

Gel Running Tray 250 X 130 mm (1 No.)

Connecting Cord (2 Nos.)

13 Well Combs 1.5 mm (1 No.)

Lid (1 No.)

Accessories

S No.	Catalog No.	Description
1	GENECO-H-MIDIP-25	Gel Running tray 25 x 13cm
2	GENECO-H-MIDIP-25UV	UV transparent Gel Running tray 25 x 13 cm (1/pack)
3	GENECO-H-LT-27	Levelling table for gel casting 20 x 27 cm
4	GENECO-H-MIDIP-EL	Spare electrode 2 Nos.

Combs

S No.	Code	Description
1	GENECO-H-MIDIP-1.5-13	Acrylic Comb MINI 13 Well 1.5mm
2	GENECO-H-MIDIP-6.0-3	Acrylic Comb MINI 3 Well 6 mm
3	GENECO-H-MIDIP-3.0-3	Acrylic Comb MINI 3 Well 3 mm
4	GENECO-H-MIDIP-CS	Comb Stand

Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com