



# GENAXY VERTICAL SYSTEMS

## MIDI Vertical Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.

### GENECO-V Vertical Systems



GENECO offers a complete range of Vertical Gel Units that combines the ease of use with high resolution separations. REGULAR as well as DUAL Models are available to run two gels simultaneously.

#### Features:

- Compatible with all type of Precast gels of major brands
- Rapid gel casting & loading
- Low buffer volumes required
- Main unit is of Acrylic & Combs & spacers are made up of Teflon

### Technical Specification

Model	GENECO-V-MIDI	GENECO-V-MIDI DUAL
Principal Material	Acrylic	Acrylic
Dimensions (L X B X H)	25 X 20 X 27 cm	23 X 18 X 27.5 cm
No. of Samples	13	26
Gel Size	16 x 14 cm	16 x 14 cm (2 Gel)
Buffer Volume	150 + 200 ml	250 + 500 ml
Combs	13 Well – Teflon 0.5 mm, 1 mm, 1.5 mm (1)	13 Well – Teflon 0.5 mm, 1 mm (2)
Teflon Spacers	0.5 mm, 1 mm, 1.5 mm (4)	0.5 mm, 1 mm (4)
Platinum Electrode	Red & Black, 1 each	Red & Black, 1 each
Connecting Cord	Red & Black, 1 each	Red & Black, 1 each
Lid	1 No.	1 No.
Glass Plates	Notch/Rectangular (2) Set	Notch/Rectangular (2) Set
Gasket	Fixed	Fixed
Clamp/Screws	1 Sets	4 Sets
Water Circulation	No	Yes

#### Accessories

S No.	Catalog No.	Description
1	GENECO-V-MIDI-GP/NP	Glass Plates (Notched or Regular)
2	GENECO-V-MIDI-Caster	Mini Gel Caster for 2 Gels
3	GENECO-V-MIDI-GDT	Gel Developing Tray with gel holder 16 x 14 cm
4	GENECO-V-MIDI-EL	Spare Electrode 2 Nos.
5	GENECO-V-MIDI-CC	Connecting cords( 2Nos.)
6	GENECO-V-MIDI-GPS	Glass Plates Stand for Midi System

#### Combs & Spacers

S No.	Code	Description
1	GENECO-V-MIDI-x-13	Teflon Combs MIDI 13 Well
2	GENECO-V-MIDI-x	Teflon Spacers (6 Nos/Pack)
3	GENECO-V-MIDI-3.0	Teflon Spacers (6 Nos/Pack)

#### Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com