

NEW unique all **GENSleek**

Soft Comfort
The smooth plunger mechanism offers least resistance for ease and effortless accuracy

Comfortable Finger Rest
For resting of the pipette on either hand

Tip Ejector
For smooth ejection of the tip

Ergonomy
Best in-hand grip for comfort

Precise Setting
Of the selected volume ensures accurate aspiration & dispensing of the liquid

Light Weight
Pipette makes your job less tiresome & improves efficiency

Ideal Shaft Shape
Helps reach contents in narrow or deep tubes

Durable Tip Cone
Offers highest resistance to shocks and chemical corrosion



Fully autoclavable UV resistance, single channel pipettes

- ❑ GENSleek is made of mechanically robust and chemically resistant materials, making it a durable and stable pipette giving consistent performance for long time.
- ❑ Single hand pipetting can be done smoothly to improve operation efficiency.
- ❑ Hand temperature doesn't permeate easily through the body of the pipette. Therefore the accuracy of the volume measurement is not affected.
- ❑ Slim nozzle can be easily inserted into tubes to aspirate, which is useful at lab work.

GENSleek, VARIABLE VOLUME PIPETTE

| CAT NO | DESCRIPTION | Volume μ l |
|-----------------|--|----------------|
| GENSleek-2 | Digital Micro Pipette, Increment 0.002 μ l | 0.1-2.5ul |
| GENSleek-10 | Digital Micro Pipette, Increment 0.02 μ l | 0.5-10ul |
| GENSleek-20 | Digital Micro Pipette, Increment 0.02 μ l | 2-20ul |
| GENSleek-50 | Digital Micro Pipette, Increment 0.1 μ l | 5-50ul |
| GENSleek-100 | Digital Micro Pipette, Increment 0.1 μ l | 10-100ul |
| GENSleek-200 | Digital Micro Pipette, Increment 0.2 μ l | 20-200ul |
| GENSleek-1000 | Digital Micro Pipette, Increment 1 μ l | 100-1000ul |
| GENSleek-5000 | Digital Micro Pipette, Increment 10 μ l | 500-5000ul |
| GENSleek-10000 | Digital Micro Pipette, Increment 20 μ l | 1000-10000ul |
| GENSleek-T-5ML | Tips for GenSleek-5000, Clear | 1-5ml |
| GENSleek-T-10ML | Tips for GenSleek-10000, Clear | 1-10ml |
| GEN-SAS-01 | Carousel Stand for 6 single channel pipettes | - |

GENSleek, FIXED VOLUME PIPETTE

| CAT NO | DESCRIPTION | Volume μ l |
|-------------------|-----------------------|----------------|
| GENSleek-FV-2.5 | Digital Micro Pipette | 2.5 μ l |
| GENSleek-FV-5 | Digital Micro Pipette | 5 μ l |
| GENSleek-FV-10 | Digital Micro Pipette | 10 μ l |
| GENSleek-FV-20 | Digital Micro Pipette | 20 μ l |
| GENSleek-FV-50 | Digital Micro Pipette | 50 μ l |
| GENSleek-FV-100 | Digital Micro Pipette | 100 μ l |
| GENSleek-FV-200 | Digital Micro Pipette | 200 μ l |
| GENSleek-FV-250 | Digital Micro Pipette | 250 μ l |
| GENSleek-FV-500 | Digital Micro Pipette | 500 μ l |
| GENSleek-FV-1000 | Digital Micro Pipette | 1000 μ l |
| GENSleek-FV-2000 | Digital Micro Pipette | 2000 μ l |
| GENSleek-FV-5000 | Digital Micro Pipette | 5000 μ l |
| GENSleek-FV-10000 | Digital Micro Pipette | 10000 μ l |