



CSQ-FC



IEF Unit

Catalog No.	Description	Price (Rs.)
<b>Accessories For CSQ33</b>		
CSQ-FC	Flexi Caster for CSQ20 & CSQ33	51,168.00
CSQ33-NG	Glass plates, Pk/2 Notched	19,866.00
CSQ33-PG	Glass plates, Pk/2	10,988.00
CSQ33-S0.25	Spacer set 0.25mm	4,588.00
CSQ33-S0.35	Spacer set 0.35mm	4,588.00
CSQ33-S1	Spacer set 1mm	4,588.00
CSQ33-S1.5	Spacer set 1.5mm	4,588.00
CSQ33-0.25-48	Comb 48 sample, 0.25mm thick, Shark tooth	6,960.00
CSQ33-0.25-96	Comb 96 sample, 0.25mm thick, Shark tooth	11,760.00
CSQ33-0.35-48	Comb 48 sample, 0.35mm thick, Shark tooth	6,960.00
CSQ33-0.35-96	Comb 96 sample, 0.35mm thick, Shark tooth	11,760.00
CSQ33-1-48	Comb 48 sample, 1mm thick, Square Tooth	6,724.00
CSQ33-1-80	Comb 80 sample, 1mm thick, Square Tooth	6,724.00
CSQ33-1.5-48	Comb 48 sample, 1.5mm thick, Square Tooth	6,724.00
CSQ33-1.5-80	Comb 80 sample, 1.5mm thick, Square Tooth	6,724.00
CSL-FHS	Fan heater sensor kit for large format vertical units CSQ20 and CSQ33	1,29,254.00
CSL-MGR	Small Glass Plate Rack	6,124.00
CSL-LGR	Large Glass Plate Rack	9,276.00

### Flatbed Isoelectric Focusing System

The redesigned CSL-IEF has been optimised to perform first-dimension isoelectric focusing (IEF) with IPG (immobilised pH gradient) strips quickly, easily and reproducibly. A high-capacity focusing tray that accommodates up to twelve IPG strips. Uniform cooling of the platform is achieved using a large ceramic plate that may be connected to any standard chiller via quick-fit connectors for enhanced strip cooling.

Catalog No.	Description	Price (Rs.)
CSL-IEF	Flatbed IEF system for IEF strips and gels	2,16,000.00
CSL-CHILLER	Chiller for Electrophoresis systems	POR
<b>Accessories</b>		
CSL-RHYDTRY	Rehydration Tray	8,880.00
CSL-FOCUSTRAY	Focusing tray with adjustable electrodes	50,400.00
CSL-IEFFRME	Replacement electrode frame	26,160.00
CSL-IEFPOS	Replacement positive electrode	29,040.00
CSL-IEFNEG	Replacement negative electrode	29,040.00
CSL-IEFPLT	Replacement Glass platform	POR
CSL-FTELECPOS	Focusing tray with adjustable electrode Positive	32,400.00
CSL-FTELECNEG	Focusing tray with adjustable electrode Negative	32,400.00

### 2-D Combination Packages

Catalog No.	Description	Price (Rs.)
CSL-IEF-KIT	1-D Combination package, includes CSL-IEF, CSL-Chiller and EV2320	POR
CSL-IEF-KIT-MINI	CSL-IEF-KIT and CVS10DSYS	POR
CSL-IEF-KIT-MINIWIDE	CSL-IEF-KIT and VS10WDSYS and IPG Converter Kit	POR
CSL-IEF-KIT-WAVE	CSL-IEF-KIT and WAVEDSYS and IPG Converter Kit	POR
CSL-IEF-KIT-MAXIPLUS	CSL-IEF-KIT and VS30DSYS	POR