





Inoculating Needle/Loops

Combination provides choice to spread or pick colonies and cultures

Transfer bacterial cultures without risk of contamination.

- Use loop end to uniformly streak a gel surface
- Use needle end for removal of colonies
- Choose between flexible polypropylene or the more

rigid ABS versions - Lot-to-lot traceability - Disposable

Item No.	Color	Material	Loop L	L x D in	L x D cm	Packaging	Sterile	UOM
HS81121A	Blue 	Flexible PP	1	7.9 x .2	20 x .4	20 ea/bag	9	20 ea/bag, 50 bags/cs 1000 ea/pk, 6 pks/cs
HS81121B		Rigid ABS						
HS81121C	Blue 	Flexible PP	10	7.9 x .2	20 x .4	20 ea/bag	9	20 ea/bag, 50 bags/cs 1000 ea/pk, 6 pks/cs
HS81121D		Rigid ABS						