

## PCR WORK STATION

### GENAXY PCR-MINI/ PCR-MIDI/PCR-MIDI-HF

▪ **PCR WORK STATION MIDI WITH HEPA FILTER (Designed for UV Irradiation of PCR Chemicals)**

<p><b>Specifications</b></p>	<ul style="list-style-type: none"> <li>• UV tubes : 1 nos.</li> <li>• Fluorescent tube : 2 nos.</li> <li>• working area LX H X D - 600mm x 500 mm x 450 mm</li> </ul> <p><b>Cabinet : MS powder coated with SS working top</b></p> <ul style="list-style-type: none"> <li>• Pipette holder : 1 No for 5 Pipettes</li> <li>• HEPA FILTER : 10 X 10 Cm 2 Nos</li> <li>• Switches at side</li> <li>• Inside Socket for running Mini Equipments</li> <li>• Sliding Door from Front side</li> </ul>
<p><b>Advantages</b></p>	<ul style="list-style-type: none"> <li>• Spacious</li> <li>• Easily transferable</li> <li>• Great design</li> <li>• Economical</li> </ul>



GEN-PCR-MIDI-F

▪ **Midi PCR Workstation (Designed for UV Irradiation of PCR Chemicals)**

<p><b>Specifications</b></p>	<ul style="list-style-type: none"> <li>• UV tubes : 1 nos.</li> <li>• Fluorescent tube : 2nos.</li> <li>• Working area LX H X D - 650 mm x 450 mm x 300 mm</li> <li>• Switches at side</li> <li>• Inside Socket for running Mini Equipments</li> <li>• Sliding Door from Front side</li> <li>• Cabinet : MS powder coated with SS working top</li> <li>• Pipette holder</li> </ul>
<p><b>Advantages</b></p>	<ul style="list-style-type: none"> <li>• Spacious</li> <li>• Easily transferable</li> <li>• Great design</li> <li>• Economical</li> </ul>



GEN-PCR-MIDI

▪ **Mini PCR Workstation (Designed for UV Irradiation of PCR Chemicals)**

<b>Specifications</b>	<ul style="list-style-type: none"> <li>• UV tubes : (1 No for mini)</li> <li>• Fluorescent tube : 1no.</li> <li>• Working area LX H X D - 250 mm x 120 mm x 120 – Mini (Cat# - GEN - PCR-MINI)</li> <li>• Switches at side</li> <li>• Sliding Door from Front side</li> <li>• Cabinet : MS powder coated with SS working top</li> </ul>
<b>Advantages</b>	<ul style="list-style-type: none"> <li>• Spacious</li> <li>• Easily transferable</li> <li>• Great design</li> <li>• Economical</li> </ul>



**GEN-PCR-MINI**

MODEL	GEN-PCR-MINI	GEN-PCR-MIDI	GEN-PCR-MIDI-HF
<b>SPECIFICATION</b>			
<b>UV tubes</b>	1 nos.		
<b>Fluorescent Tube</b>	1 nos.	2 nos.	
<b>Working Area LX H X D (MM)</b>	250mm x 120 mm x 120 mm	650mm x 450 mm x 300 mm	600mm x 500 mm x 450 mm
<b>Cabinet</b>	MS Powder Coated with SS Working Top		
<b>Pipette Holder</b>	---	---	1 No for 5 Pipettes
<b>HEPA FILTER</b>	---	---	10 X 10 Cm 2 Nos.
<b>Switches</b>	At Side		
<b>Inside Socket</b>	---	For Running Mini Equipments	
<b>Sliding Door</b>	From Front side		

**ORDERING INFORMATION**

CATALOG NO.	DESCRIPTION	QTY
GEN-PCR-MINI	PCR WORK STATION MINI	EACH
GEN-PCR-MIDI	PCR WORK STATION MIDI	EACH
GEN-PCR-MIDI-HF	PCR WORK STATION MIDI WITH HEPA FILTER	EACH

**Genaxy Scientific Pvt. Ltd.**