

GENAXY Pipette Controller

This pipette fillers excellent ergonomic design, the precise control of pipetting speed & Large LCD display provides visual confirmation of remaining battery charge and speed settings. It allows your pipetting work to take it leisurely and impressively, easy to add and transfer the culture medium during cell seeding and cell passage, gives you complete pipetting comfort and speed up the experimental process.

Features

- Large LCD display provides visual confirmation of remaining battery charge and speed settings Light weight and ergonomic design
- Pipetting speed can be controlled simply and adjusted steplessly.
- Adjustable motor speed (3 part motor speed for choice)
- Pipette Fillers for use with pipettes from 0.1-100mL
- Efficient Lithium polymer rechargeable battery ensures number of dispensing is approx. 2,000 (with a 25 mL pipette)
- without recharging



Catalog No.	Description	Volume
GENPET-ePF	Pipette Controller with Large LCD Display, Double safety valve and Hydrophobic PTFE Membrane Filter 0.45µm	1-100 ml

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com