

# PLASMID DNA ISOLATION KITS

Kits for isolation of high quality Plasmid DNA ; complete range of standard Mini Prep, Midi Prep, Maxi Prep kits and 96-well Kits. High recovery of plasmid DNA : Yield from mini is 30ug , from Midi 750ug and from Maxi kit 1.5mg. These kits are time saving, versatile and safe. They eliminates the use of phenol, chloroform, and CsCl and minimizes the exposure to hazardous chemicals.


## Plasmid Mini, Midi and Maxi Kit

Catalog No.	Description	Kit Size
FAPDE 100	FavorPrep™ Plasmid DNA extraction mini kit	100 Preps
FAPDE 300	FavorPrep™ Plasmid DNA extraction mini kit	300 Preps
FAPDE 002	FavorPrep™ Plasmid DNA extraction midi kit	25 Preps
FAPDE 002-1	FavorPrep™ Plasmid DNA extraction midi kit	50 Preps
FAPDE 002-EF	FavorPrep™ Endotoxin Plasmid DNA extraction midi kit	25 Preps
FAPDE 003	FavorPrep™ Plasmid DNA extraction Maxi kit	10 Preps
FAPDE 003-1	FavorPrep™ Plasmid DNA extraction Maxi kit	20 Preps
FAPDE 003-EF	FavorPrep™ Endotoxin Plasmid DNA extraction Maxi kit	10 Preps



## FavorPrep™ 96 Well Plasmid Kit

Catalog No.	Description	Kit Size
FAPWE 96004	FavorPrep™ 96 Well Plasmid DNA extraction Kit for High throughput	96 X 4 Plates



## Plasmid Extraction Kit for QIAcube

Catalog No.	Description	Kit Size
FAPDE-AUTO-024	Plasmid DNA Extraction Mini Kit for QIAcube Automatic Machine	24 Preps

## FavorFilter Plasmid DNA Extraction Maxi Kit

- Sample size: Up to 240ml of Bacterial Culture.
- Expected Yield: up to 500ug DNA.
- Purified DNA can be used for Transfection, In Vitro Transcription and Restricted Enzyme digestion.

Catalog No.	Description	Kit Size
FAFTE 001	FavorFilter Plasmid DNA Extraction Maxi Kit	4 Preps
FAFTE 001-1	FavorFilter Plasmid DNA Extraction Maxi Kit	10 Preps
FAFTE 001-EF	FavorFilter Endotoxin Free Plasmid DNA Extraction Maxi Kit	4 Preps
FAFTE 001-1-EF	FavorFilter Endotoxin Free Plasmid DNA Extraction Maxi Kit	10 Preps

