

# Gel Casters—Small Format Vertical Electrophoresis

Enables several gradient or homogenous gels to be cast simultaneously



## Multiple Gel Casters

- Pour uniform gels from the top; pump gradient gels into an inlet port at the bottom of the unit
- Sandwiches fill at the same time at the same rate from the same solution. The result: identical gels.
- Glass plate sandwiches seal leak free without the use of messy grease or tape
- Once cast, gels can be wrapped in plastic and refrigerated for several weeks

## Ordering Information

Cat. #	Description
SE215	Ten Gel Caster for 10 x 8 cm plates Includes: 10 x 8 cm Rectangular Glass Plates—20 pcs 10 x 8 cm Notched Alumina Plates—10 pcs Spring Clamps—2 pcs
SE275	Four Gel Caster for 10 x 8 cm plates Includes: 10 x 8 cm Rectangular Glass Plates—10 pcs 10 x 8 cm Notched Alumina Plates—4 pcs Spring Clamps—2 pcs
SE235	Four Gel Caster for 10 x 10.5 cm plates Includes: 10 x 10.5 cm Rectangular Glass Plates—5 pcs 10 x 10.5 cm Notched Alumina Plates—4 pcs Spring Clamps—4 pcs

### All three above include:

- Caster Body with Gasket
- Face Plate with Inlet Port
- Rectangular Glass Plates—size and qty varies
- Notched Alumina Plates—size and qty varies
- Spring Clamps—qty varies
- Wax Paper—100 sheets
- Space-Saver Plate
- Polycarbonate Filler Sheets—5 pcs
- Black Vinyl Cap
- Filler Plug Set

Order SE211A combs and spacers separately, see page 31.

## Dual Gel Caster

- Cast one or two gels
- Accommodates both 10 x 8 cm and 10 x 10.5 cm plates
- Uses the same rugged and effective clamp/cam casting method as used with the SE600 and SE400 Series casters

## Ordering Information

Cat. #	Description
SE245	Dual Gel Caster

### Includes:

- Casting Cradle
- Casting Clamp Assemblies—2 pcs
- Sealing Gasket Set
- Cams—4 pcs

Order glass plates, SE211A combs, and spacers separately, see pages 22 and 31.

