

UV Transilluminators

Available in single and dual wavelength formats, in 21x21cm and 21x26cm sizes. Standard features include a high/low intensity safety switch and an efficient starter that allows each of the six 8-Watt UV tubes to energize quickly without flickering, while special filter glass minimizes unwanted background light.

Catalog No.	Description	Price (Rs.)
CSLUVTS254	UV Transilluminator, small, 21 x 21 cm, 254 nm	2,04,000.00
CSLUVTS312	UV Transilluminator, small, 21 x 21 cm, 312 nm	2,04,000.00
CSLUVTS365	UV Transilluminator, small, 21 x 21 cm, 365 nm	2,04,000.00
CSLUVTSDU0	UV Transilluminator, small, 21 x 21 cm, 254/365 nm	2,18,400.00
CSLUVTSDU0312	UV Transilluminator, small, 21 x 21 cm, 254/312 nm	2,18,400.00
CSLUVTSDU0365	UV Transilluminator, small, 21 x 21 cm, 312/365 nm	2,18,400.00
CSLUVTS254L	UV Transilluminator, large, 21 x 26 cm, 254 nm	2,47,200.00
CSLUVTS312L	UV Transilluminator, large, 21 x 26 cm, 312 nm	2,47,200.00
CSLUVTS365L	UV Transilluminator, large, 21 x 26 cm, 365 nm	2,47,200.00
CSLUVTSDUOL	UV Transilluminator, large, 21 x 26 cm, 254/365 nm	2,16,630.00
CSLUVTSDU0312L	UV Transilluminator, large, 21 x 26 cm, 254/312 nm	2,16,630.00

Accessories

CSL-T254	8W, Replacement UV bulb 254nm	8,100.00
CSL-T365	8W, Replacement UV bulb 365nm	8,100.00
CSL-UVPS22	UV Transpatent Cutting Platform 22 x 22 cm	11,340.00
CSL-UVPS27	UV Transpatent Cutting Platform 22 x 27 cm	12,310.00
CSL-F21254	Replacement Filter 21cm - 254nm	1,20,000.00
CSL-F21365	Replacement Filter 21cm - 365nm	1,20,000.00
CSL-F26254	Replacement Filter 26cm - 254nm	1,58,400.00
CSL-F26312	Replacement Filter 26cm - 312nm	1,58,400.00
CSL-F26365	Replacement Filter 26cm - 365nm	1,58,400.00
CSL-C	Replacement transilluminator cover	16,560.00
CSL-UVSCRN	UV to white light Transilluminator screen converter	73,720.00

Blue Light Transilluminator


Cleaver's new SAFEBlue LED transilluminator offers the user a safe way to view and document their samples. SAFE : No DNA damage to samples, No UV Light. Using High Purity LED light. Fast start up. Strong Metal enclosure with stainless steel filter frame.

Catalog No.	Description	Price (Rs.)
SAFEVIEW	Blue light transilluminator, 21x21cm	POR
SAFEVIEW-MINI2	BLUE LIGHT TRANSILLUMINATOR 15.3 x 15.3cm	92,350.00
proBLUEVIEW	proBLUEVIEW mini transilluminator, 12x18cm, with flat pack hood	POR
RVMINI	Mini Blue Light transilluminator, 11.2x7.4cm with filter	POR

duoVIEW Transilluminator


Dual UV & Blue light technology in one transilluminator. Available in Single or dual wavelength models. Industrial leading long life UV filter. Fast start up.

Catalog No.	Description	Price (Rs.)
DUOVIEW312	UV 312 & BLUE Light transilluminator, 21x21 cm	POR
DUOVIEW254	UV 254 & BLUE Light transilluminator, 21x21 cm	POR
DUOVIEW365	UV 365 & BLUE Light transilluminator, 21x21 cm	POR
DUOVIEW254/312	UV 254/312 & BLUE Light transilluminator, 21x21 cm	POR
DUOVIEW254/365	UV 254/365 & BLUE Light transilluminator, 21x21 cm	POR



UV Transilluminator

Filter Size	21x21cm 21x26cm
Light Source	8Wx6 TUBES
UV Resistant Plastic cover	33x25 cm (WxD)
Unit dimensions (WxLxH)	34x29.5x10cm



Blue Light Transilluminator

Filter Sizes	21 x 21cm
Light Source	470nm Blue LED's
Size	340x270x130mm
Weight	5kg
Voltage	110-240V (Selectable)



duoVIEW Transilluminator

Filter Sizes	21 x 21cm
Light Source	470nm Blue LED's
OR	Single Wave 8W x 5 tubes Dual Wave 8W x 9 tubes
Intensity Switch	High (100%)/Low (70%) single
Size	400x190x350mm
Weight	10.5kg