


Cell Culture Dishes

1HVW

Features

- IP
-
-
-
-
-
-



Tissue Culture Dishes For Cell Attachment

Cat.No.	Spec(mm)			Cell Growth Area (cm ²)	Recommended Medium Volume (mL)	/Pack	/Case
	Type	Height	Diameter				
706001	35	12	33.6	8.5	1.8-2.7	20	500
New 706201	35 _A (With Gripping Ring)	12	33.6	8.5	1.8~2.7	20	500
705001	60	15	55.2	22.9	4.2-6.3	20	500
704004	100	20	87.2	57.6	11-16.5	5	300
704002	100	20	87.2	57.6	11-16.5	10	300
704001	100	20	87.2	57.6	11-16.5	20	300
704202	100 _A (With Gripping Ring)	20	86.9	59.3	11-16.5	10	300
704201	100 _A (With Gripping Ring)	20	86.9	59.3	11-16.5	20	300
715001	150	25	139.9	150.1	30.4-45.6	5	100