

CENTRIFUGE

Product Features

- Electronic locks. Centrifugation only when the lid is locked.
- Capacitive touchpad, responsive.
- LCD display. Real-time display of speed, time and operating status.
- The lid opens automatically after centrifugation.
- Centrifugal time and speed adjustable. Adjustable speed during centrifugation.
- Suitable for 0.2mL / 0.5mL / 1.5mL / 5mL centrifuge tubes without change the rotor.
- Fast Lifting Speed.
- Low Noise



Rotors



Universal Centrifuge



5 mL Centrifuge



High - Speed Centrifuge

Technical Specifications

Catalog No.	Name	Range Of Rotation	Max RCF	Rotor	RPM & TIME	SIZE	Voltage
C10609	Universal Centrifuge	500-6000 RPM	2400g	General Rotor (0.2/0.5/1.5/5.0ml)	Adjustable	260x220x150	110-220V 50-60HZ
C10605	5ml Centrifuge	500-6000 RPM	2400g	Rotor for 5mL centrifuge tubes	Adjustable	260x220x150	110-220V 50-60HZ
C11202	high - speed centrifuge	500-12000 RPM	8800g	Rotor for 1.5mL & 2.0mL centrifuge tubes	Adjustable	260x220x150	110-220V 50-60HZ