



● Internal Thread Cryogenic Vials

- Used for storage in the gas phase of liquid nitrogen
- Special screw cap designed with silicone O-ring for a secure seal
- Many sizes available as self-standing with universal locking base

Cat. No.	Volume(ml)	Description	Sterile	Qty./Pack	Qty./Case
618101	0.5	Internal, round bottom, self-standing	+	50	500
606101	1.2	Internal, round bottom, self-standing	+	50	500
607101	2.0	Internal, round bottom, self-standing	+	50	500
608101	4.0	Internal, round bottom, self-standing	+	50	1,000
609101	5.0	Internal, round bottom, self-standing	+	50	1,000



● Nest Colored Inserts

- 6 colored cap inserts for easy identification

Cat. No.	Color	Qty./Pack	Qty./Case
611001	White	100	1,000
611002	Red	100	1,000
611003	Yellow	100	1,000
611004	Blue	100	1,000
611005	Green	100	1,000
611006	Purple	100	1,000