

Elisa Plates (Immunoplates)



Product Features

- Designed for ELISA method
- Medical grade polystyrene
- Conform to SBS/ ANSI
- Both detachable and undetectable plates are available
- Flat and transparent well bottom
- High binding (protein binding: 300-400 ng/cm²)
- Medium binding (protein binding: 200-300 ng/cm²)


Undetachable ELISA Plate

- Uniform thickness and stripe well size, flat

bottom , CV <5%

Detachable ELISA Plate

- Clear alphanumeric codes for better test flexibility and identification
- Special designed stripe-holder for eliminating the change of disorientation
- Designed for easy handling and automation
- 1*8 stripes design for easy readouts
- Excellent well uniformity, CV <5%



Cat.No.	Description	Spec (well)	Color	/Pack	/Case
504201	Detachable, High Binding, Clear	96	White bottom & Transparent well	5	50
514201	Undetachable, High Binding, Clear	96	Transparent	5	50