

Mini Vortex Mixer

NEST Mini Vortex Mixer is a test tube oscillator that focuses on the jog function. The trigger-type start switch is safe and convenient. It is one of the indispensable tools in the laboratory. It can be used to oscillate test tubes and centrifuge tubes.

Product Features

- ❑ Extremely compact and easy operation;
- ❑ Powerful vortexing for tubes up to 30mm;
- ❑ Head is made of durable material to ensure long life;
- ❑ Refined steel base is designed to ensure product stability.



Technical Specifications

Product Name	Mini Vortex Mixer
Cat. No.	105003
Oscillation Mode	Circumference Oscillation
Orbital Diameter	Φ4.5mm
Max. Capacity	80mL
Max. Tube Diameter	Φ30mm
Speed	3000rpm
IP rating (DIN EN 60529)	IP40
Ambient Temp.	5~40°C
Outside Dimensions	Φ96xH85mm
Power	5W
Power Supply	AC100~240V,50/60Hz
Weight	0.55kg