

customer satisfaction through business ethics....

PETG Erlenmeyer Flasks



- Feature
- PETG Bottle & HDPE Cap
 - Sterilized by E-Beam, Nonpyrogenic
 - Individually wrapped in sterile bags
 - Vent caps with 0.22um hydrophobic filters for gas exchange without contamination.
 - Non-toxic and environmentally friendly.

Cat.No.	Volume (mL)	Cap Style	Size(mm)			Qty./Case
			Height	Bottleneck Diameter	Bottom Diameter	
781001	125	Plug Seal	106.8	34	66	24
781011	125	Vent	106.8	34	66	24
782001	250	Plug Seal	137.8	34	83	12
782011	250	Vent	137.8	34	83	12
783001	500	Plug Seal	175.8	39	101	12
783011	500	Vent	175.8	39	101	12
784001	1000	Plug Seal	213.3	39	127	6
784011	1000	Vent	213.3	39	127	6

Marketed By : **GENAXY SCIENTIFIC PRIVATE LIMITED**, Janakpuri, New Delhi 110058;
Phones : +91-11-4761 9999 (100 Lines) Email: technical@genaxy.com

