

## BETS-M1 Orbital Decoloring Shaker

Super engineering ABS plastic molding that is elegant and made for long-lasting. It is also environmental-friendly and made without any spray paint.

The plastic platform is resistant from acid base chemicals and easy to wash. The plastic will not change shape and does not leak. The unique design of the on/off button, the speed knob, the timing knob, and the power cord is fenced with plastic that is resistant from liquids, therefore avoiding any kind of short-circuit power off.

The digital display and the timing device is easy to read. The mini-design is a space saver and can be stored in the freezer.



### Technical Specification

Model	BETS-M1
Frequency	40 – 240 rpm
Time Range	0 – 120 min
Rotation Radius	15 mm
Display	Digital
Capacity	30 W
Power	220 V
Plate Dimensions	240 x 230 mm
Overall Dimensions	240 x 290 x 130 mm
Weight	5 kg