

KB Series

The KB series intelligent shakers are an economic and environmentally friendly compact shakers for many laboratory applications. The units have small footprints suitable for stable operation on the laboratory bench. Digital display and setting of speed and time by push button benefits convenient handling. Different variants of KB series allow rocking, shaking or 3D motion for routine laboratory application like gel staining, washing etc.



KB-700



KB-800



KB-900

Technical Specification

Model	KB-700	KB-800	KB-900
Power	220V	220V	220V
Capacity	35W	35W	35W
Frequency	0 - 80 turns / min	0 - 80 rpm	20 - 240 rpm
Speed control	Digital	Digital	-
Rocking range	-	Up - down 20 mm	-
Rotation range	-	-	Turning radius of 15 mm
Vibrating type	3 - D motion	-	-
Timing range	0 - 99.99 hours	-	-
Plate	320x265mm	280x260mm	320x265mm
Overall dimensions	370x335x190mm	350x260x150mm	370x335x145mm