

QB-8002 Micro-Plate Mini Shaker

New design and elegant, this mini shaker is made specifically for two 96 well micro plates.

It is also made for long - lasting.

The motor will run `continually without overheating.

ABS plastic molding that is classic and made to last .

It is also environmental friendly and made without any spray paint.

Various voltage available and works between 80 - 240 Voltage.

The platform is designed to fit one or two micro plates.



Technical Specification

Model	QB-8002
Power	80/240V 50/60HZ
Capacity	20 W
Rotation radius	15 mm
Frequency	0 - 1500 rpm
Speed	Step-less speed regulation
Oscillation radius	2 mm
Platform type	QB-8002/2micro-plate
Overall Dimensions	260 x 150 x 110 mm
Weight	3 kg