

QB - 328 Ultra - Small Rotation Incubator

Ultra - small design with a simple and beautiful structure.

Suitable for use in low temperature cabinet and constant temperature cabinet.

Angle of rotation can be adjusted and the horizontal position can be used for small mixing while the vertical position can be used for complete mixing.

Rotation type speed regulation can arbitrarily set the mixing requirements of 10 - 80 rev / min.

You can choose to buy three kinds of tube fixtures with different sizes, which are suitable for different sample cells. It is used to process the suspended biology or the blood sample which requires for gas exchange.

It is equipped with simple fixtures which can be quickly installed or removed.



Technical Specification

Model	QB - 328
Power supply	220V
Capacity	40 W
Temperature range	4°C-60°C
Mixing angle	Can be adjusted freely
Rotational velocity	5 - 80 turns / min
Rotate speed control	Knob type with step - less speed regulation
Overall dimensions	240×210×310mm
Standard configuration	1.5ml*48 or 15ml*20 or 50ml*12 choose one of three