

Vortex-6 Vortex Shaker

Shaking or Vortexing Action: Variable speed control allows slow speed shaking action up to high-speed vortexing. In addition, each attachment is designed to produce different mixing actions, allowing users to choose the correct action for each application.

Does not walk: Tough metal housing combined with “just enough weight”, provides a stable platform for all types of mixing.

Unmatched Reliability: Many years of proven performance in laboratories all over the world.

Timer Function: timer is available

Rugged and reliable, most popular vortex mixer. Variable speed control allows high speed shaking or vortexing with Speed (RPM)2800 ,Product Weight :4.8 Kg Dimensions (base) :170 x 120 x 170mm.



Technical Specification

Model	Vortex-6
Power	220V
Capacity	50W
Rotation Speed	2800 rpm
Control	Light - operated induction
Induction Distance	30cm
Working	Continuous touch and speed regulation
Plate	Exchangeable between bowl type and flat type
Overall Dimensions	170×120×170mm
Weight	4.8KG