

ROTATOR

AG / Adjustable speed and angle of rotation

- Durable construction for long-term operation

HAG / Use with a wide variety of accessories


- Optional rocking platform for washing

Common facts



- Various speed and adjustable mixing angle 0°~90°
- Continuous or timed operation
- Wide choice of accessories to be used, including conical tubes, bottles and 250ml flasks
- Add up to 2 disks together for application
- LED display (shown actual speed or time)

► An examples to use

- Bacteria / Yeast culture (small scale)
- Genomic DNA preparation
- Plasmid preparation



► HAG



► AG



Specification

Model	AG	HAG
Angle Range	0° ~ 90°	
Speed Range	7 ~ 70rpm	
Timer (speed)	99hour 59min	
Temp.Resist	Lower than 65°C(Dry air)	
Weight	7.5kg	8kg
Dimensions (mm)	W220xD236xH246	W510xD236xH246
Wattage	35W	
Power	Free Voltage (AC100~240V)	
Frequency	50/60Hz	

Disk

Model Number	Description
D50-8	Disk / 50ml conical tube-capacity : 8ea
D15-16	Disk / 15ml conical tube-capacity : 16ea
D1.5-32	Disk / 1.5ml centrifuge tube-capacity : 32ea
D1.5-60	Disk / 1.5ml centrifuge tube-capacity : 60ea
D250-4	Disk / 250ml Flask-capacity: 4ea
PO-4	Post for disks


Roller Drum

Model Number	Description
RD16-151	Roller Drum / Up to 15mm Diameter tube-capacity : 151ea
RD18-104	Roller Drum / Up to 17mm Diameter tube-capacity : 104ea
RD25-42	Roller Drum / Up to 29mm Diameter tube-capacity : 42ea

Rotisseries (for HAG)

Model Number	Description
12RT-4	80mm Diameter 200mm length bottle-capacity : 4ea
12RT-6	60mm Diameter 200mm length bottle-capacity : 6ea
12RT-10	40mm Diameter 200mm length bottle-capacity : 10ea
12RT-8V	50ml Conical tube-capacity : 8ea-vertical mount
12RT-16V	50ml Conical tube-capacity : 16ea-vertical mount
12RT-24H	50ml Conical tube-capacity : 24ea-horizontal mount
12RT-200-64	1.5ml Centrifuge tube-capacity : 64ea
12RT-200-96	1.5ml Centrifuge tube-capacity : 96ea
12RT-200-128	1.5ml Centrifuge tube-capacity : 128ea
12RT-300-120	1.5ml Centrifuge tube-capacity : 120ea
12RT-300-180	1.5ml Centrifuge tube-capacity : 180ea
12RT-300-240	1.5ml Centrifuge tube-capacity : 240ea
12RK-HAG	Rocking Platform(W285xD195xH65, tilt angle : ±7°)

ACCESSORIES (AG / HAG)



D1.5-32

D1.5-60

D15-16

D50-8

D250-4

RD16-151

RD18-104

RD25-42

PO-4

ACCESSORIES (HAG)



12RT-4

12RT-6

12RT-10

12RT-8V

12RT-16V

12RT-24H

12RT200-64

12RT-200-96

12RT200-128

12RT300-120

12RT300-180

12RT300-240

Silicon cap Bottle

12RK-HAG