

orbital / reciprocating SHAKER

► SH30



■ SH30 / Use with a wide variety of large capacity platforms

- Smooth and powerful orbital motion at low speed

■ SH30L / Large capacity platforms - smooth reciprocating motion

- Smooth and powerful reciprocating motion at low speed


► Common facts

- Continuous or timed operation
- Optional second platform available to double the working area, large upper platform for a variety of flasks and tubes
- Other type flask platform also available upon special request
- Standard elastic ties for easy fitting of flasks and conical tubes
- Non-slip mat as standard, optional dimpled pad prevents tubes from rolling around

► SH30L



- An examples to use
- Bacteria / Yeast culture



Specification

Model	SH30	SH30L
Orbital range	14mm	-
Reciprocate range	-	30mm
Speed range	40 ~ 300 rpm	
Motor	DC Motor	
Timer (speed)	180min	
Temp. Resist	Lower than 65°C(Dry air)	
Loadable weight	8kg	
Weight	15kg	
Dimensions (mm)	W330xD330xH157	W330xD330xH162
Platform (mm)	W300xD300	
Wattage	35W	
Power	Free Voltage (AC100-240V)	
Frequency	50/60Hz	

Accessories

Model Number	Description
DP-300	Dimpled Pad (W300xD300mm)
TY300	Second platform (with posts : 4ea)
UP400	400x400 mm Upper Platform
UP500	500x500 mm Upper Platform
SR300	Flexible silicon wire rack (W300xD300mm)

Other Accessories

Model Number	Description
CP50	50ml flask holder
CP100	100ml flask holder
CP250	250ml flask holder
CP500	500ml flask holder
FT302	Universal Table for flask holder
FT302-50-13	Flask Table for 50ml flask-Capacity : 13ea
FT302-100-13	Flask Table for 100ml flask-Capacity : 13ea
FT302-250-8	Flask Table for 250ml flask-Capacity : 8ea

ACCESSORIES

