



GENTip Dual-layered Filters

GENTip Dual-layered Filters

GENTip dual-layered filter has utilized the superb hydrophobicity of UHMW-PE to prevent aerosols from passing through filter and entering into pipette, also protect against impurities of pipette from contaminating samples.

Maximum Air-permeability, Strengthen Blocking Capacity

GENTip dual-layered filter consists of two layers with different pore size. The pore size of the upper blue layer is $5\mu\text{m}$, and the pore size of the down white layer is $10\mu\text{m}$, which assures the consistent air flow and strengthen the blocking capacity of aerosols.

Dual-layered Filters, Dual Barriers

With dual barriers : the $10\mu\text{m}$ white layer blocking droplets, splashes and suspending particles; the $5\mu\text{m}$ blue layer preventing biomolecules efficiently, the dual-layered filter has reinforced blocking property.

GENTip dual-layered filter is free of DNase/RNase, endotoxins and PCR inhibitors for liquid handling including PCR, radioactive, biotoxic, chemical corrosive and volatile samples.

10 μl Dual Aerosol-Barrier Pipet Tips

Catalog	Volume	Color	Packaging	Type	Qty./Pack	Price/Pack	Price/10Pack (Rs.)
GEN-DFB-10	0.5-10 μl	Natural	Bulk	Non-Sterilized	1000 tips	3,395.00	32,255.00
GEN-DFB-10-S	0.5-10 μl	Natural	Bulk	Pre-Sterilized	1000 tips	4,085.00	38,815.00
GEN-DFB-10XL	0.5-10 μl	Natural	Bulk	Non-Sterilized	1000 tips	3,535.00	33,585.00
GEN-DFB-10XL-S	0.5-10 μl	Natural	Bulk	Pre-Sterilized	1000 tips	4,245.00	40,325.00
GEN-DFB-10-R-S	0.5-10 μl	Natural	Racked	Pre-Sterilized	10 racks	5,465.00	51,915.00
GEN-DFB-10XL-R-S	0.5-10 μl	Natural	Racked	Pre-Sterilized	10 racks	5,605.00	53,245.00

200 μl Dual Aerosol-Barrier Pipet Tips

Catalog	Volume	Color	Packaging	Type	Qty./Pack	Price/Pack	Price/10Pack (Rs.)
GEN-DFB-200	1-200 μl	Natural	Bulk	Non-Sterilized	1000 tips	3,395.00	32,255.00
GEN-DFB-200-S	1-200 μl	Natural	Bulk	Pre-Sterilized	1000 tips	4,085.00	38,815.00
GEN-DFB-200-R-S	1-200 μl	Natural	Racked	Pre-Sterilized	10 racks	5,465.00	51,915.00



1000 μl Dual Aerosol-Barrier Pipet Tips

Catalog	Volume	Color	Packaging	Type	Qty./Pack	Price/Pack	Price/10Pack (Rs.)
GEN-DFB-1000	100-1000 μl	Natural	Bulk	Non-Sterilized	1000 tips	4,345.00	41,275.00
GEN-DFB-1000-S	100-1000 μl	Natural	Bulk	Pre-Sterilized	1000 tips	5,210.00	49,495.00
GEN-DFB-1000-R-S	100-1000 μl	Natural	Racked	Pre-Sterilized	6 racks	3,775.00	35,865.00

