

**GENAXY<sup>®</sup>**  
Customer Satisfaction Through Business Ethics



**GENAXY Mini-6K & 10K Centrifuge**

# GENAXY Mini-6K & 10K Centrifuge

## Features


- Compact shape, colorful, streamlined design, small size, mute technology, ultra-low noise operation
- Suitable for a variety of centrifuge tubes, one machine for multiple purposes
- Compliant with international safety regulations.
- Safety design: flip-open switch function, the cover will automatically stop when it is opened.
- The outer cover is made of flame-retardant composite material, which is not easy to break
- Micro DC motor and carbon fiber rotor
- This new Mini Centrifuges accommodates 6 x 1.5-2.0ml tubes & 2 x 8 x 0.2ml strip tubes.
- These units are ideal for quick spin downs, microfiltration & cell separations available in 6,000 rpm and 10,000rpm.
- Two Models to choose from
- Minifuge for 1.5-2.0ml tubes and 8 x 0.2ml strip tubes (available in 6K model only)

GEN-Mini-6K is supplied with 3 rotors:

- 6 x 1.5/2.0ml angle rotor
- 2 x 8 x 0.2ml strip rotor
- Rotor for slide

GEN-Mini-10K is supplied with 1 rotor:

- 6 x 1.5/2.0ml angle rotor



## Technical Specification

Model	GEN-MINI-6K	GEN-MINI-10K
Power requirement	30W	40W
Relative centrifugal force (RCF)	2,000 x g	5,000 x g
Max. rotational speed	6,000 rpm 6x1.5/2.0ml angle rotor	10,000 rpm
Rotors	2 - 8x0.2ml strip rotor Rotor for slide	6x1.5/2.0ml angle rotor
Dimensions	175 x 148 x 118 mm	175 x 148 x 118 mm
Weight	0.7kg	0.7kg

## Accessories

Code	Description
AS-08061-01	Allen wrench
AS-08061-02	1.5/2.0ml rotor
AS-08061-03	Strip tube rotor
AS-08061-04	Slide rotor
AS-08061-05	Adapters for 0.5-0.6ml tubes
AS-08061-06	Adapters for 0.2ml tubes

