


Septa Mat



GENAXY Septa Mat is designed to be used on 3500 series genetic analyzer. Made of research-grade silicone, they can be autoclaved, and following a bleach wash and ethanol rinse protocol, they can also be reused.

- These sealing mats minimize evaporation in a wide range of temperatures in PCR and storage applications
- Made of research-grade silicon
- Offers unsurpassed protection from Volatile solvents
- Silicon pre-slit well caps permit a pipette tip or probe to enter the well without damage to the cap

Cat. Number	Description	Qty. / Case
GEN-96-SEPTA-3100	96-well Silicon Septa Mat for 3500 series genetic analyzer	50

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com