

## PCR plate sealing mat



Sealing mats minimize evaporation in a wide range of temperatures in PCR and storage applications. Made of research-grade silicone, they can be autoclaved, and following a bleach wash and ethanol rinse protocol, they can also be reused.

- 96 well and 384-well PCR Sealing Mat
- For use with 96 round well PCR microplates Like:  
GEN-9-010-PCR-NS-C, GEN-9-020-PCR-NS-C,  
GEN-9-010-PCR-FS-WC, GEN-9-020-PCR-NSE-C,  
GEN-9-020-PCR-SS-C, GEN-9-020-ES-ABC, GEN-9-010-ES-ABC
- For use with 384 round well PCR microplates Like:  
GEN-3-PCR-FS-C, GEN-3-00-FS-C-ABI

| Cat. Number     | Description                               | Qty. / Case           |
|-----------------|-------------------------------------------|-----------------------|
| GEN-PCOMAT-96W  | Sealing Mat for 96-well PCR Plate, White  | 10 pcs/pk ; 5 pk/case |
| GEN-PCOMAT-384W | Sealing Mat for 384-well PCR Plate, White | 10 pcs/pk ; 5 pk/case |

**Genaxy Scientific Pvt. Ltd.**