


Discovering excellence with Innovation

DNA Ladders

The DNA ladder is ready-to-use, requiring no additional preparation or dilution, making it convenient and time-saving. It is mixed with tracking dye for easy visualization during gel electrophoresis, ensuring monitoring. Offering reliable sizing, it provides precise and reproducible DNA fragment analysis, enabling effective molecular weight determination. This ladder is suitable for agarose gel electrophoresis and other gel systems, making it a versatile and dependable tool for a variety of laboratory applications



Ordering Information

Ready to Use: DNA ladders are fully ready to load on gel and also contains tracing dye

Sharp and Intense Bands: with defined DNA amount in each band its sharp and intense

Applications: for agarose gel electrophoresis and DNA -PAGE

CATALOG NO	Inclusions	Pack size
GBI-100201	100bp DNA ladder	500µl
GBI-100202	100bp Plus DNA ladder	500µl
GBI-100208	100bp extended DNA ladder	500µl
GBI-100203	200bp DNA Ladder	500µl
GBI-100205	250bp plus DNA ladder	500µl
GBI-100206	500bp DNA ladder	500µl
GBI-100207	1kb DNA ladder	500µl
GBI-100208	1kb plus DNA ladder	500µl

More Information

✉ info@genaxy.com

🌐 www.genaxy.com

📍 GENAXY SCIENTIFIC PVT LTD
Near Janak Cinema, Janakpuri,
New Delhi