

CELL CULTURE PLATES



GENAXY cell culture plates have round edges for better handling. Low profile, low-evaporation lid design assures uniform growth. Choose tissue culture-treated or non-treated surfaces. Non-treated growth surface is more hydrophobic than tissue culture surface and is similar to the surface of a bacteriological Petri dish. Small feet on the bottom corners of each plate eliminate scraping of the microscope lens.

- Sterilized by E-beam
- High clarity, 100% virgin polystyrene
- Vacuum Plasma tissue culture treatment
- Clear lot number for batch traceability
- Peel-to-open, individually wrapped
- Markings of well coordinates available for 12, 24, 48, 96 and 384 well plates
- Non-Pyrogenic, DNase/RNase FREE

Cat. Number	Description	Recommended Volume	Cell Growth Area	Qty. / Pack
GEN-TCP-6S	6 Well with LID, Sterile	3	9.5	50
GEN-TCP-12S	12 Well with LID, Sterile	2	3.6	50
GEN-TCP-24S	24 Well with LID, Sterile	1	1.9	50
GEN-TCP-48S	48 Well with LID, Sterile	0.5	0.95	50
GEN-TCP-96S	96 Well with LID, Sterile	0.1-0.2	0.32	100

Genaxy Scientific Pvt. Ltd.