


CELL CULTURE FLASKS



GENAXY cell culture flasks are available in cell growth areas ranging from 25cm² to 225cm². Tissue culture treated or non-treated flasks available with a vent cap or plug seal cap to meet your needs

- Packaged in sterile, zip-sealable bags
- High clarity, 100% virgin polystyrene
- Frosted writing and clear graduations
- Notched bottom for slip free stacking
- Sterilized by E-Beam, Vent caps (with breathable .22um membrane)
- Clear lot number for batch traceability
- Non-Pyrogenic, DNase/RNase FREE

Cat. Number	Description	Recommended Volume	Cell Growth Area	Qty. / Pack
GEN-TCF-25-PS	Flask, 50ml, Sterile, Plug cap	5-7.5	25	200
GEN-TCF-25-VS	Flask, 50ml, Sterile, Vent cap	5-7.5	25	200
GEN-TCF-75-PS	Flask, 250ml, Sterile, Plug cap	15-22.5	75	100
GEN-TCF-75-VS	Flask, 250ml, Sterile, Vent cap	15-22.5	75	100
GEN-TCF-175-PS	Flask, 750ml, Sterile, Plug cap	35-52.5	175	40
GEN-TCF-175-VS	Flask, 750ml, Sterile, Vent cap	35-52.5	175	40
GEN-TCF-225-PS	Flask, 950ml, Sterile, Plug cap	45-67.5	225	25
GEN-TCF-225-VS	Flask, 950ml, Sterile, Vent cap	45-67.5	225	25

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com