

CELL CULTURE DISHES



GENAXY cell culture dishes are ideal for all types of culturing where microscopic examination is required. Numeric indicators on the bottom of each plate allow users to identify the location of cells.

- Ultra pure virgin polystyrene, optically clear conforming to USP class IV
- Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment for excellent cell adherence, optimum and uniform growth of adherent cells
- Vented ports for effective gaseous exchange
- Easy to handle and better gripping
- Stable notched bottom design prevents spillage and contamination.
- Stackable for easy storage and handling
- Sterilized by gamma irradiation Non-cytotoxic
- Non-Pyrogenic, DNase/RNase and Human DNA free

Cat. Number	Description	Recommended Volume	Cell Growth Area	Qty. / Pack
GEN-TCD-35-S	35mm Dishes, with LID, Sterile	1.8-2.7	9	500
GEN-TCD-60-S	60mm Dishes, with LID, Sterile	4.2-6.30	21	500
GEN-TCD-100-S	100mm Dishes, with LID, Sterile	11-16.50	55	300
GEN-TCD-150-S	150mm Dishes, with LID, Sterile	30.4-45.6	152	100

Genaxy Scientific Pvt. Ltd.