



# PETRI DISHES

*TRANSPARENCY  
MATTERS*

Genaxy petri dishes are made of Medical Grade virgin polystyrene and is best for microbiological purpose. These plates comes with or without ventilation ribs. It also reduces the condensation during incubation. Flared lid skirt and squared corners allows easy one hand operation. Various packaging options are available in the form of single wrap, double wrap, sleeve of 10 and sleeve of 20.


## Unique Features

- » Genaxy petri dishes are molded from medical grade virgin polystyrene
- » Manufactured under ISO Class 8 clean room facility
- » USFDA 21 CFR Compliant
- » Clear top & bottom Viewing surfaces
- » Radiation sterilized : Gamma or ETO
- » Certificates of Sterilization processing available on request
- » Available in Various packaging: individual wrap, sleeve of 10 pieces, & sleeve of 20 pieces. Sleeves can be single or double wrap.
- » Certified DNase/RNase, Human DNA, and Pyrogen Free
- » Different types of petri plate are available to serve different applications:
- » **90mm vented petri dishes** comes with Ventilation ribs for free air circulation & reduces condensation during incubation
- » **90mm non-vented dishes** have flat lid that prevents evaporation , encouraging long incubation period and anaerobic applications. Consistent flatness allows for even media distribution
- » **Two & Three compartment petri dishes** for easy identification which is specifically designed for simultaneous studies of organism in a common environment.
- » **Media saver petri dishes** comes in smaller diameter thus uses lesser media. Ideal for use in automation. Available in 35mm & 60 mm diameter
- » **Contact/RHODAC plates** are useful for bio-burden counts on working surface. Comes with Domed base or flat base with grid molded on the surface of the base. Compatible with air samplers.

*Contact plates are now available with "LOCKING SYSTEM"*



# Genaxy Petri Dishes



Flared lid skirt for easy one hand operation especially with gloves

Raised edges on lid for easy stacking

Ventilation ribs allow for free air circulation and reduces condensation during incubation

Consistent flatness allow for even media distribution



FOR MORE PRODUCT INFORMATION  
SCAN THE QR CODE



## Petri Dishes

Cat. Number	Description	Qty. / Pk
GEN-PTD-35	35x12mm Petri Dishes, PS, Individual Peel Off, "MEDIA SAVER", Gamma Sterile	500
GEN-PTD-35-10	35x12mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, "MEDIA SAVER", Gamma Sterile	500
GEN-PTD-60	60x15mm Petri Dishes, PS, Individual Peel Off, "MEDIA SAVER", Gamma Sterile	300
GEN-PTD-60-10	60x15mm Petri Dishes, "MEDIA SAVER", PS, Sleeve of 10, Gamma Sterile	300
GEN-PTD-90-DL	90x15mm Petri Dishes, PS, Individual Peel Off, Fully stackable lid, Gamma Sterile with Locking Feature	460
GEN-PTD-90-DL-10	90x15mm Petri Dishes, PS, Packed in sleeve of 10, Fully stackable lid, Gamma Sterile with Locking Feature	540
GEN-PTD-90	90x15mm Petri Dishes, PS, Individual Peel Off, Fully stackable lid, Gamma Sterile	460
GEN-PTD-90PP	90x15mm Petri Dishes, Individual Peel Off, Fully stackable lid, Polyprop-lene, Autoclavable	200
GEN-PTD-90-ETO	90x15mm Petri Dishes, PS, Individual Peel Off, ETO Sterile	460
GEN-PTD-90-10	90x15mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, Gamma Sterile	540
GEN-PTD-90-10-NV	90x15mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, Gamma Sterile, Non-Vented	540
GEN-PTD-90-10-ETO	90x15mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, ETO Sterile	540
GEN-PTD-90-10-ETO-NV	90x15mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, ETO Sterile, Non-Vented	540
GEN-PTD-90-10-D/W	90x15mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, Double Wrapped, Gamma Sterile	540
GEN-PTDA-90-10	90mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, Automation Compatible	540
GEN-PTD-2C-90	90x15mm Petri Dishes, PS, 2 Compartment, Individual Peel Off, Gamma Sterile	460
GEN-PTD-3C-90	90x15mm Petri Dishes, PS, 3 Compartment, Individual Peel Off, Gamma Sterile	460
GEN-CTP-65	65x15mm Contact Dish, Convex Base with Grid	400
GEN-CTP-65-D	65x15mm Contact Dish, Convex Dome Base with Grid	400
GEN-CTP-65-DL	65x15mm Contact Dish with lock, Convex Dome Base with Grid	400
GEN-PTD-100-20	100x20mm Petri Dishes, PS, Packed in sleeve of 10 Pcs, Gamma Sterile	300
GEN-PTD-150	150x15mm Petri Dishes, PS, Individually Packed, Gamma Sterile	100
GEN-PTD-130SQ	130x130x15mm Petri Dishes, PS, Square, Gamma Sterile	150

## Genaxy Scientific Pvt. Ltd.