

EXTERNAL THREAD CRYO VIALS



- These Cryogenic vials are manufactured from medical grade polypropylene
 - Withstand temperatures to -196 Degree C.
 - All vials can be color-coded with cap inserts, which are sold separately.
 - Vials have graduation marks and a large marking spot for easy labeling of contents.
 - Caps have O-Ring for leakproof sealing.
 - Self-standing bottom Design.
 - All tubes are Gamma sterile.
 - Shatter proof and pressure resistant vials.
 - ISO 9001:2015, 13485:2016 & 14001:2015 certified.
- DNase and RNase - free and non-pyrogenic.
 - Do not use cryo vials in the liquid phase of liquid nitrogen.

Cat. Number	Description	Qty. / Pack
GEN-CVE-1.2ml	1.2ml Cryogenic Vial, Graduated, Sterile.	500
GEN-CVE-2.0ml	2.0ml Cryogenic Vial, Graduated, Sterile.	500
GEN-CVE-5.0ml	5.0ml Cryogenic Vial, Graduated, Sterile.	250
GEN-CVE-1.2ml-X	1.2ml Cryogenic Vial, Amber Color, Graduated, Sterile.	500
GEN-CVE-2.0ml-X	2.0ml Cryogenic Vial, Amber Color, Graduated, Sterile.	500
GEN-CVE-5.0ml-X	5.0ml Cryogenic Vial, Amber Color, Graduated, Sterile.	250

Genaxy Scientific Pvt. Ltd.

INTERNAL THREAD CRYO VIALS

Cat. Number	Description	Qty. / Pack
GEN-CVI-1.2ml	1.2ml Cryogenic Vial, Graduated, Sterile.	500
GEN-CVI-2.0ml	2.0ml Cryogenic Vial, Graduated, Sterile.	500
GEN-CVI-5.0ml	5.0ml Cryogenic Vial, Graduated, Sterile.	250
GEN-CVI-1.2ml-X	1.2ml Cryogenic Vial, Amber Color, Graduated, Sterile.	500
GEN-CVI-2.0ml-X	2.0ml Cryogenic Vial, Amber Color, Graduated, Sterile.	500
GEN-CVE-5.0ml-X	5.0ml Cryogenic Vial, Amber Color, Graduated, Sterile.	250

CAP INSERTS FOR INTERNAL THREAD CRYO VIALS

Cat. Number	Description	Qty. / Pack
GEN-BC-CI-Y	Cap Inserts for Internal Thread Vials, Yellow.	50
GEN-BC-CI-L	Cap Inserts for Internal Thread Vials, Lime Green.	50
GEN-BC-CI-P	Cap Inserts for Internal Thread Vials, Pink.	50
GEN-BC-CI-PU	Cap Inserts for Internal Thread Vials, Pink.	50

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com