

# PASTEUR PIPETTES

—TRANSFER WITH EASE—




Genaxy pasteur pipettes are excellent for use in all general laboratory applications, including blood banks, hematology, chemistry, and bacteriology. Molded from translucent, low density Polyethylene, which is inert to biological fluids and most acids. These are unbreakable, Non-toxic, all-in-one pipettes which eliminate the hazard of broken glass and exposure to infectious material. It can be sealed and refrigerated. Low-binding surface reduces the loss of cells and valuable proteins.

## Unique Features

- » Manufactured under ISO Class 8 Clean room facility
- » Molded from translucent, low density Polyethylene, which is inert to biological fluids and most acids
- » Unbreakable, Non-toxic, all-in-one pipettes eliminate the hazard of broken glass and exposure to infectious material
- » Low-binding surface reduces the loss of cells and valuable proteins
- » Well-suited for blood banking.
- » Available Volume: 0.2ml, 1ml, 3ml, 5ml and 10ml
- » Sterile or Non-Sterile both options are available
- » Sterile products are ETO sterile
- » DNase/RNase and Pyrogen free
- » US FDA 21 CFR compliant



# Genaxy Pasteur Pipettes



FOR MORE PRODUCT  
INFORMATION  
SCAN THE QR CODE



## Pasteur Pipettes

Cat. Number	Description	Qty. / Pack
GEN-151	1ml Pasteur Pipette, Non-Sterile	1000
GEN-151-S	1ml Pasteur Pipette, Individually Packed, Sterile	100 X 5
GEN-153	3ml Pasteur Pipette, Non-Sterile	1000
GEN-153-S	3ml Pasteur Pipette, Individually Packed, Sterile	100 X 5
GEN-155	5ml Pasteur Pipette, Non-Sterile	1000
GEN-155-S	5ml Pasteur Pipette, Individually Packed, Sterile	100 X 5
GEN-150	0.2ml Pasteur Pipette, Non-Sterile.	1000
GEN-156	10ml Pasteur Pipette, Non-Sterile.	500
GEN-156-S	10ml Pasteur Pipette, Individually Packed, Sterile.	500

**Genaxy Scientific Pvt. Ltd.**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058  
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : [info@genaxy.com](mailto:info@genaxy.com)